

WEBINAR

# HEAVY METALS TESTING: Considerations for an ICP Analysis of the Big Four

THURSDAY, NOVEMBER 11  
9:00–10:00AM EST

SPEAKERS:

Autumn Phillips  
R&D Chemist



Lesley Owens, PhD  
Technical Support  
Manager



# Heavy Metals Testing: Considerations for an ICP Analysis of the Big Four

- Overview of Heavy Metals Testing
  - Impacted Industries
  - Current Regulations
  - Current Events
  - Future
- Overview of Issues
  - Instrumental
  - Working solution(s)
  - Sample Prep
- Summary
  - IV Products
  - Technical Support



# Overview of Heavy Metal Testing

- Impacted Industries
  - Everyone!
- Current Regulations
  - USP 232
  - US EPA
  - Consumer Product Safety Act (CPSA)
    - ASTM Standard F963-17 [<25 ppm As, <60 ppm Hg, <75 ppm Cd, and <90 ppm Pb]
  - FDA
    - Hg in fish [1 ppm] & cosmetics [1 ppm]
    - Color additives in cosmetics [3 ppm As]



# Overview of Heavy Metal Testing

- Current Regulations (continued)
  - Mercury-Containing and Rechargeable Battery Management Act
    - Button Cell batteries [25 ppm Hg]
    - Rechargeables [Pb & Cd]
  - RoHS
    - EU directive that is shaping policy in the US
    - 10 substances regulated [Cd < 100 ppm; Pb & Hg < 1000 ppm]
  - Model Toxics in Packaging Legislation
    - Focuses on product packaging [Cd, Hg, & Pb < 100 ppm by weight]
  - CA Proposition 65
    - List of restricted chemicals
    - Maximum concentrations vary by product





# Overview of Heavy Metal Testing

- Current Events
  - Baby rice cereal
    - Voluntary recall
    - Samples exceed limit for inorganic arsenic (100 ppb)
    - October 8, 2021
  - Furniture
    - Recall
    - Surface paint exceeds the limit for Pb (90 ppm)
    - October 28, 2021
- Future – more testing!



# Overview of Issues

- Instrumental – spectral and mass interference(s); washout
- Working solution(s) – stability and compatibility
- Sample Preparation




# Instrumental Issues: Arsenic

- Has no cationic chemistry (As(III) will exist as arsenite, As(V) will exist as arsenate)
- The oxidation state of As **DOES** influence the intensities obtained using ICP-OES.
- There are inconsistencies regarding which form of As, (III) or (V), gives higher intensities on ICP-OES, but in our observations, As(III) gives higher intensities.



# Instrumental Issues: Arsenic




- The green line is 10ppm As(III)
- The purple line is 10ppm As(V)



# Instrumental Issues: Arsenic

- It is important to note the starting material used for your As standards ( $\text{As}$  metal,  $\text{As}_2\text{O}_3$ ,  $\text{As}_2\text{O}_5$ ), as this will affect the form of As in solution.
- Ensure that the As species is the same in your standards as in your samples ( $\text{As(III)}$  or  $\text{As(V)}$ )



# Instrumental Issues: Arsenic

Atomic Spectroscopic Information: (red text indicates severe at ~ concs.)

Technique / Line	Estimated D.L.*	Order	Type	Interferences
ICP-OES 189.042 nm	0.05/.005 µg/mL	1	atom	Cr
ICP-OES 193.696 nm	0.1/.01 µg/mL	1	atom	V, Ge
ICP-OES 228.812 nm	0.1/.01 µg/mL	1	atom	Cd, Pt, Ir, Co
ICP-MS 75 amu	30 ppt	n/a	M <sup>+</sup>	<sup>40</sup> Ar <sup>35</sup> Cl, <sup>59</sup> Co <sup>16</sup> O, <sup>36</sup> Ar <sup>38</sup> Ar <sup>1</sup> H, <sup>38</sup> Ar <sup>37</sup> Cl, <sup>36</sup> Ar <sup>39</sup> K, <sup>150</sup> Nd <sup>2+</sup> , <sup>150</sup> Sm <sup>2+</sup>

\*ICP-OES D.L.'s are given as radial / axial view

- The 189.042nm wavelength is usually the best line.
- Always check multiple lines to make sure they agree.
- Major Cd interference on 228.812nm wavelength. Do not use this line if Cd is present.
- If chloride is present, there will be a major interference from ArCl on ICP-MS
  - Use collision cell technology (He mode), if possible, to break up chloride interferences.



# Instrumental Issues: Mercury

- **Hg MUST be in the  $\text{Hg}^{+2}$  form for accurate ICP measurements!**
- Hg sticks to plastic in  $\text{HNO}_3$  matrices → many intro system components are composed of plastic!
- Glass intro system preferred for Hg
- Initial analyses should not be affected, however memory effects due to Hg washing out over time can affect subsequent analyses.
  - Thiourea or higher HCl rinses can help with washout issues
- Best stabilized in HCl. Chloride will cause interferences for As and Se on ICP-MS
  - Use collision cell (He mode) to break up interferences



# Instrumental Issues: Mercury

- $\text{Hg}_2^{2+}$  can be converted to a mix of  $\text{Hg}^{2+}$  and  $\text{Hg}^0$  in the presence of chloride or other ligands.
- If Hg is present, even in small part, as the metallic ( $\text{Hg}^0$ ) species, the nebulization/transport efficiency will be significantly higher for Hg and signal could increase several hundred percent.
- Extremely high recoveries → reduced Hg species
- Low results → Hg loss due to adsorption or precipitation of  $\text{Hg}_2^{2+}$  with chloride





# Instrumental Issues: Mercury

Atomic Spectroscopic Information: (red text indicates severe at ~ concs.)

Technique / Line	Estimated D.L.*	Order	Type	Interferences
ICP-OES 184.950 nm	0.03/.005 µg/mL	1	atom	
ICP-OES 194.227 nm	0.03/.005 µg/mL	1	ion	V
ICP-OES 253.652 nm	0.1/.03 µg/mL	1	atom	Ta, Co, Th, Rh, Fe, U
ICP-MS 202 amu	9 ppt	n/a	M <sup>+</sup>	<sup>186</sup> W <sup>16</sup> O

\*ICP-OES D.L.'s are given as radial / axial view

## Additional preferred lines:

ICP-OES: 435.835 nm

ICP-MS: 199amu, 200amu, \*201amu (best line if W is in the sample, do not use if Re is in sample)



# Instrumental Issues: Cadmium

Atomic Spectroscopic Information: (red text indicates severe at ~ concs.)

Technique / Line	Estimated D.L.*	Order	Type	Interferences
ICP-OES 214.438 nm	0.003/.0003 µg/mL	1	ion	Pt, Ir
ICP-OES 228.802 nm	0.003/.0003 µg/mL	1	atom	Co, Ir, <b>As, Pt</b>
ICP-OES 226.502 nm	0.003/.0003 µg/mL	1	ion	Ir
ICP-MS 111 amu	11 ppt	n/a	M <sup>+</sup>	<sup>95</sup> Mo <sup>16</sup> O

\*ICP-OES D.L.'s are given as radial / axial view

\*Major As interference on 228.802nm wavelength. Do not use this line if As is present.

- Other recommended masses: 112 and 114amu



# Instrumental Issues: Lead

Atomic Spectroscopic Information: (red text indicates severe at ~ concs.)

Technique / Line	Estimated D.L.*	Order	Type	Interferences
ICP-OES 168.215 nm	0.03/0.003 µg/mL	1	ion	Co
ICP-OES 220.353 nm	0.04/0.006 µg/mL	1	ion	Bi, Nb
ICP-OES 217.000 nm	0.09/0.03 µg/mL	1	atom	W, Ir, Hf, Sb, Th
ICP-MS 208 amu	5 ppt	n/a	M*	<sup>192</sup> Pt <sup>16</sup> O, <sup>192</sup> Os <sup>16</sup> O

\*ICP-OES D.L.'s are given as radial / axial view

- Additional masses: 204 (**Hg**, Os, and Yb interference), 206 (Pt and Os interferences), and 207amu (Ir interference)



# Instrumental Issues

- Salting out may be an issue for high salt matrices.
  - Cone conditioning prior to analysis may be necessary
  - Or you may need to dilute samples further
  - Aerosol dilution can be used, if available, to dilute sample before introducing to the plasma
- Corrosion of instrument components due to high acid
- Clogging of the nebulizer
  - Contact instrument manufacturer regarding TDS limitations and other limitations of your particular instrument





# Instrumental Issues

- **Matrix matching** samples and standards is very important for achieving accurate results.
  - Account for spectral interferences
  - Account for signal suppression from high matrix
  - Ensure consistent nebulization efficiency between samples and standards
- If the solutions and standards can not be matrix-matched, **standard additions** is a great method to account for matrix effects.
  - More information: <https://www.inorganicventures.com/icp-guide/standard-addition-internal-standardization-and-isotope-dilution>



# Instrumental Issues

- **Internal Standards** can help correct nebulization or plasma related effects. Choose an internal standard that:
  - Is not contained in the sample (lanthanides are common, but should be avoided if fluoride is present)
  - Is stable in the solution (compatible with matrix and other elements)
  - Gives sufficient signal to noise ratio
  - Behaves similarly in the plasma
  - Will not have interferences from other elements in the sample or interfere on other analytes
  - Has a mass similar to the analyte of interest (for ICP-MS)
    - If there are multiple analytes, multiple internal standards may be necessary to cover the desired mass range.
    - The most common internal standards are  $^6\text{Li}$ , Sc, Y, In, Tb and Bi.
- <https://www.inorganicventures.com/catalogsearch/result/?q=internal%20standard>




# Instrumental Issues

- **Reference Materials:** use a well-characterized material that is similar to your sample that can be used to validate your method
  - Are you getting the correct results for each analyte in the RM using your method?
  - Is anything being lost?
- **Run Method Blanks!** This will help determine whether contamination is being introduced.



# Instrumental Issues: Washout

- Of the Big Four, **Hg** is going to be the one that has washout issues
- For OES, we recommend reconditioning the spray chamber surface with RBS™-25 (2.5%) → will make the surface less “sticky”
  - Not recommended for MS (high Na)
- Sample tubing is the main issue for Hg sticking (especially PVC tubing)
- HCl / Thiourea
  - 1-10% v/v HCl
  - 0.5% w/v Thiourea
  - Also good for Au
- NH<sub>4</sub>OH
  - 1-5% v/v for OES or MS
- \*If Triethanolamine (TEA) has been used for HF neutralization in a prior analysis, make sure to switch out your waste solution before Hg analysis as Hg can be reduced to Hg<sup>0</sup> and travel back up the lines to the instrument causing memory effects





# Working Solutions Issues

- Pb is a common contaminant in pipet tips and acid reagents ( $\text{HNO}_3$ )
- Keep in **acid** matrix!
- Avoid **neutral/basic** media
  - Many cationic metals may form insoluble arsenates when mixed with As under pH neutral conditions.
  - Cd can form insoluble carbonate and hydroxide.
  - Pb can form insoluble carbonate, borate, sulfate, sulfite, sulfide, phosphate, oxalate, chromate, tannate, iodate, and cyanide
  - Hg can form insoluble carbonate
- Avoid mixing Pb with  **$\text{Cr}^{+6}$**  and  **$\text{H}_2\text{SO}_4$**
- Avoid mixing Hg with **tartrate**, as Hg can be reduced to the metallic form.
- Cadmium chloride, bromide, and iodide are soluble in water.
  - \* $\text{CdI}_2$  is one of the few iodides soluble in ethanol
- **All** Cd compounds are soluble in excess NaI, due to the formation of the complex ion,  $\text{CdI}_4^{2-}$ .



# Working Solutions Issues

- As and Pb, along with Bi, can have issues when mixed at higher concentrations.
- If a solution contains any two of the following elements, As, Bi, or Pb, at concentrations greater than 500ppm, a higher acid matrix of at least 10% v/v HNO<sub>3</sub> may be needed to maintain stability.
- If As and Bi are both present at concentrations of 1000ppm or greater, then at least 15% v/v HNO<sub>3</sub> is required.
- However, high acid concentrations can cause Pb to drop out as the nitrate, so be careful with the order of additions and do not add Pb to a concentrated HNO<sub>3</sub> environment.
  - Same issues exist for Ba and Cs.
- Pb can also precipitate in HCl. Pb at 500ppm or greater may need higher HCl for stability (recommend 20% v/v HCl or greater).



# Working Solutions Issues

- Reducing environments can convert  $\text{Hg}^{2+}$  to the  $\text{Hg}_2^{2+}$  dimer or the metallic form.
- $\text{Hg}_2^{2+}$  disproportionation into metallic mercury and  $\text{Hg}^{2+}$ :  $\text{Hg}_2^{2+} \xrightleftharpoons{\text{Cl}^-} \text{Hg}^0 + \text{Hg}^{2+}$
- $\text{Hg}^{2+}$  can reduce to  $\text{Hg}_2^{2+}$  if  $\text{NO}_2$  gases remain in the container from reaction/production process
- Adding  $\text{HNO}_3$  and boiling can convert  $\text{Hg}_2^{2+}$  back to the  $\text{Hg}^{2+}$  (recertification may be required)



→ Dilute with  $\text{H}_2\text{O}$  to make stock products/concentrates



- \*Be Aware that boiling a solution containing As(III) with  $\text{HNO}_3$  will oxidize As(III) to As(V)

# Working Solutions Issues



- Hg in HNO<sub>3</sub> matrix can adsorb to plastic container walls (~1ppm adsorbed)
- More of a concern for lower Hg concentrations (plastic not recommended for <200ppm Hg)
- Solution: May need to use borosilicate glass containers
- Solution concerns:
- Cannot have HF in standards contained in glass; HF is caustic to glass
- Glass has higher levels of contaminants than LDPE/HDPE
- Elements of concern: Al, Ba, B, Ca, Ga, Fe, Ni, K, Na, Sr, Zn, Zr



# Working Solutions Issues

- Hg can be stabilized with  $\text{AuCl}_3$  (generally 1ppm  $\text{AuCl}_3$ . See EPA bulletin: [https://www.inorganicventures.com/pub/media/wysiwyg/files/mercury\\_preservation\\_techniques.pdf](https://www.inorganicventures.com/pub/media/wysiwyg/files/mercury_preservation_techniques.pdf)) or HCl.
  - Be aware of chloride interference for As and Se on ICP-MS.
  - $\text{AuCl}_3$  should be included in blanks, standards, and samples.
- HCl matrix is recommended
- If Ag is present, excess chloride will need to be added to keep Ag in solution and the solution will be photosensitive.
- If Tl is present, it should be in the  $\text{Tl}^{+3}$  state as  $\text{Tl}^{+1}$  may precipitate as the chloride

\*There are also methods for stabilizing Hg with L-cysteine in dilute  $\text{HNO}_3$ , but we have not used this methods in-house yet. (See [Fresquez MR, Pappas RS, Watson CH. Establishment of Toxic Metal Reference Range in Tobacco from U.S. Cigarettes. J. Anal. Toxicol. 2013;37:298-304.](#))



# Sample Preparation: Organic Sample Containing the Big Four

- Sample preparation technique will depend on the type of sample and the analytes of interest. Do you care about everything, or can some leftover insoluble be filtered out?
  - HF may be required to dissolve everything in the sample. (ex: Si, Ti, Sn, Sb, Zr, W, Ge, Ta, Hf, Nb)
- We recommend using a closed vessel or reflux condenser to prevent loss of volatile species
- Microwave digestion technique would be best
  - Use of  $\text{HNO}_3$  and a small amount of HCl is recommended if analyzing for only As, Cd, Pb, and Hg. HCl needed to prevent Hg loss.
  - A small amount of  $\text{H}_2\text{O}_2$  may need to be added to complete oxidation of certain organic species. For analysis of these four elements, it is not necessary/recommended.
  - \*If  $\text{H}_2\text{O}_2$  is used, always add  $\text{H}_2\text{O}_2$  **LAST**, in very small increments, and proceed ***slowly*** and ***cautiously***!





# Sample Preparation: Organic Sample Containing the Big Four


- Be aware of all contamination sources (use high purity DI water and acids, clean bottles/containers (we recommend LDPE), may need to leach materials, clean workspace, materials used to grind or blend the sample, etc.)
  - Prepare method blanks
- Avoid using glass or metal containers/materials at any point in the sample preparations process (contaminants)





# Sample Preparation: Organic Sample Containing the Big Four

- Small sample size (~0.5g for cannabis)
  - Ensure that your samples are as homogenous as possible
- For sticky samples, may need to freeze before grinding
- Some possible methods:
  - EPA Methods: 200.7, 200.8, 6020A, 3052
  - USP 232/233 and ICH Q3D
  - FDA Methods 4.4 and 4.7
  - AOAC Methods 2015.01 and SMPR® 2020.001
  - Standard Methods 3030F
  - Many More!



# Sample Preparation: Additional techniques for different sample types by element

## Arsenic

- Metal: soluble in 1:1 H<sub>2</sub>O / HNO<sub>3</sub>
- Oxides: exist in crystalline and amorphous forms. The amorphous form is more water-soluble. Most oxides typically dissolve in dilute acidic solutions when boiled.
- Minerals: 1 g of powdered sample is fused in a Ni<sup>0</sup> crucible with 10 g of a 1:1 K<sub>2</sub>CO<sub>3</sub> / KNO<sub>3</sub> mix, and the melt extracted with hot water
- Organic Matrices: 0.2-0.5 g of the sample are fused with 15 g of a 1:1 Na<sub>2</sub>CO<sub>3</sub> / Na<sub>2</sub>O<sub>2</sub> mix in a Ni<sup>0</sup> crucible. The fuseate is extracted with water and acidified with HNO<sub>3</sub>.
  - Ensure that the mixture is covered with a layer of the 1:1 Na<sub>2</sub>CO<sub>3</sub> / Na<sub>2</sub>O<sub>2</sub> mixture to prevent loss of organoarsenic compounds
- As reacts with Pt, so a platinum crucible should be avoided for ashing.
- Zr crucibles should also be avoided to prevent loss of As due to formation of zirconium arsenate.
- If samples are made from arsenic trioxide, solution must be basic pH to dissolve into solution
- If standards are made from As metal, the solution must be heated with HNO<sub>3</sub> for at least 1 day after all metal has dissolved to convert all As to As<sup>+5</sup>.
  - Take extra precaution when using this method to avoid build-up of H<sub>2</sub>O<sub>2</sub> in the mixture → can lead to an explosion!

## Cadmium

- Metal: soluble in HNO<sub>3</sub>
- Oxides: soluble in HCl or HNO<sub>3</sub>
- Ores: dissolve in HCl / HNO<sub>3</sub>, then heat to fumes with H<sub>2</sub>SO<sub>4</sub>. The silica and lead sulfate are filtered off after the addition of water.
- Organic: dry ash at 450°C and dissolve ash in HCl
  - Cadmium chloride is volatile at temps over 400°C
- Alternative option: sulfuric / peroxide acid digestion:
  - Heat w/ conc H<sub>2</sub>SO<sub>4</sub> to fumes, then add 30% H<sub>2</sub>O<sub>2</sub> **dropwise** to finish oxidation
  - Make sure there is always an EXCESS of H<sub>2</sub>SO<sub>4</sub>
  - Be **patient** and proceed **slowly**!



# Sample Preparation: Additional techniques for different sample types by element

## Mercury

- Use a closed vessel or condenser for acid digestions
- Avoid ashing! (Volatility and Toxicity of  $\text{Hg}^0$ )
- Metal or oxide: dissolve in  $\text{HNO}_3$
- Organics: Extremely toxic! Please use extreme caution when working with organomercury compounds
  - Must be kept in solution with oxidizing agent to prevent loss as the metallic form!
  - Hg in biological material: heat with  $\text{H}_2\text{SO}_4$  or  $\text{HNO}_3$  or both and potassium permanganate (in excess)
    - Oxidizing agent: permanganate
    - $\text{H}_2\text{SO}_4$  alone may not fully free all Hg
    - Use a closed vessel!
    - Low heat ( $\sim 50\text{-}600^\circ\text{C}$ )
  - \*Can also use perchloric acid instead of permanganate (For tips and safe use of perchloric acid and other acid digestions, please see Section 12 of our Trace Analysis Guide: <https://www.inorganicventures.com/trace-analysis-guide/acid-digestions-of-organic-samples> )

## Lead


- Metal, Ores and Alloys: best dissolved in 1:1  $\text{H}_2\text{O}$  /  $\text{HNO}_3$
- Oxides: soluble in  $\text{HNO}_3$ , with the exception of  $\text{PbO}_2$ , which is soluble in  $\text{HCl}$  or  $\text{HF}$
- Organic Matrices: dry ash and dissolve in dilute  $\text{HCl}$ .

Do not heat when dissolving to avoid precipitation of  $\text{SiO}_2$ .



# IV Products for Heavy Metal Testing

- Stock
  - IV-STOCK-65
    - USP <232> / ICH Q3D Class 1 Oral Elemental Impurities
- Custom
  - Over 160 unique products containing As, Hg, Cd, and Pb




30 µg/mL Hg  
15 µg/mL As  
5 µg/mL Cd & Pb



# Technical Support – Available to Everyone

## Online Resources at [inorganicventures.com](http://inorganicventures.com)



The image shows a standard periodic table with elements color-coded by groups. Calcium (Ca) is highlighted in blue. A callout box for Calcium provides the following information:

- Calcium
- Atomic Weight: 40.078
- Oxidation State: +2

Below the main table is a separate row for the Lanthanides and Actinides, with Promethium (Pm) highlighted in blue.

Customers can visit our website's Tech Center, which includes:

- Interactive Periodic Table
- Sample Preparation Guide
- Trace Analysis Guide
- ICP Operations Guide
- Expert Advice
- And much, much more.

