

2022 INORGANIC VENTURES WEBINAR SERIES

How It's Made

Custom CRM Edition

Thursday | July 21 | 9:00 a.m. EST

ASHLEY JONES
TECH SUPPORT CHEMIST



Intro to Inorganic Ventures' Certified Reference Materials



ISO Accredited



ACCREDITED
ISO 17034
0883.02 BY A2LA



ACCREDITED
ISO 17025
0883.01 BY A2LA



REGISTERED
ISO 9001
QSR-1034 BY QSR



NIST traceable - We use NIST Standard Reference Materials (SRM™ and lot number specified on the certificate) for calibration and direct comparison.



Over 63,000 unique solutions available for order, 10-14 new solutions added daily



inorganicventures.com

How IV's CRMs are made

1. New custom products developed by Technical Support Chemists in a two step, PhD reviewed process.
2. Certified in-house concentrates manufactured from purity tested starting material.
 - Single or multielement high concentration CRMs used for aliquoting into our custom CRM products.



How IV's CRMs are made continued – manufacturing overview

3. Once ordered, solutions go to the manufacturing department for production

- Gravimetric preparation
- CRMs are produced under our scope of accreditation, consistency maintained through audited Work Instructions

1. Technicians obtain their manufacturing breakdowns for solutions ordered.
2. Technicians collect in-house concentrates of appropriate concentrations and verify product, concentration, lot number, and density against breakdown.
3. Aliquots are weighed out according to final concentrations requested and added to the blend – once all aliquots are added the solution is diluted to final ordered volume by weight.
4. Solutions are mixed at least 50x by inversion before retains are taken for Quality Control and solution is passed off to be packaged.



How IV's CRMs are made continued – certificate of analysis

4. Submitted to QC for verification of certified values calculated from the gravimetric manufacturing process.

- COA generated by Product Documentation department

Certificate includes:

- ISO information
- Product Description
- Cert. Values/info on assay(s) used
- NIST Traceability info
- TMI if applicable
- Additional info on storage and handling



Scan me to navigate to CoA/SDS search tool!



[inorganicventures.com](https://www.inorganicventures.com)

1.0 ACCREDITATION / REGISTRATION

INORGANIC VENTURES is accredited to ISO 17034, "General Requirements for the Competence of Reference Material Producers" and ISO/IEC 17025, "General Requirements for the Competence of Testing and Calibration Laboratories". Inorganic Ventures is also an ISO 9001 registered manufacturer (QSR Certificate Number QSR-1034).



2.0 PRODUCT DESCRIPTION

Product Code: Multi Analyte Custom Grade Solution

Catalog Number:

Lot Number:

Matrix: 5% (v/v) HNO₃
1.1% (v/v) HF

Value / Analyte(s): 100 µg/mL ea: Tin,
Molybdenum,
Tungsten,
50 µg/mL ea: Zirconium,
25 µg/mL ea: Antimony, Silver

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ANALYTE	CERTIFIED VALUE	ANALYTE	CERTIFIED VALUE
Antimony, Sb	µg/mL	Molybdenum, Mo	µg/mL
Silver, Ag	1 µg/mL	Tin, Sn	µg/mL
Tungsten, W	5 µg/mL	Zirconium, Zr	µg/mL

Density 1.028 g/mL (measured at 20 ± 4 °C)

Assay Information:

ANALYTE	METHOD	NIST SRM#	SRM LOT#
Ag	ICP Assay	3151	160729
Ag	Volhard	999c	999c
Ag	Calculated		See Sec. 4.2
Mo	ICP Assay	3134	130418
Mo	Calculated		See Sec. 4.2
Sb	ICP Assay	3102a	140911
Sn	ICP Assay	3161a	140917
Sn	Calculated		See Sec. 4.2
W	ICP Assay	3163	140606
W	Calculated		See Sec. 4.2
Zr	ICP Assay	3169	130920
Zr	Calculated		See Sec. 4.2

Benefits of gravimetric preparation

- More reliable and consistent than volumetric preparation
 - Avoid pitfalls of plastic consumables
 - adhesion
 - Increased accuracy
 - Limit affect of atmospheric pressure and temperature
 - Small margin of error (aliquots measured to 4 decimal places below 30g, 2 decimal places for +30g)
 - Viscosity
 - Ability to adjust aliquots weighed in secondary container
 - Reduce cost
 - Consumables – transfer pipettes vs pipette tips
 - Calibration of pipettes



Benefits of gravimetric preparation

- More reliable and consistent than volumetric preparation
 - Avoid pitfalls of plastic consumables
 - adhesion
 - Increased accuracy
 - Limit affect of atmospheric pressure and temperature
 - Small margin of error (aliquots measured to 4 decimal places below 30g, 2 decimal places for +30g)
 - Viscosity
 - Ability to adjust aliquots weighed in secondary container
 - Reduce cost
 - Consumables – transfer pipettes vs pipette tips
 - Calibration of pipettes

Gravimetric Vs Volumetric Aliquot

Let's make a 100mL solution of 1000ppm Ca.
(Ca concentrate has concentration of 119,910ppm & density of 1.341 g/mL).

$$\left(\frac{100\text{mL} \times 1000\text{ppm}}{119910\text{ppm}} \right) = \mathbf{0.83396\text{mL}}$$

$$0.83396\text{mL} \times \frac{1.341\text{g}}{\text{mL}} = \mathbf{1.1183\text{g}}$$

Our volumetric aliquot is 0.8340mL – this should weigh 1.1183g



Benefits of gravimetric preparation

- More reliable and consistent than volumetric preparation
 - Avoid pitfalls of plastic consumables
 - adhesion
 - Increased accuracy
 - Limit affect of atmospheric pressure and temperature
 - Small margin of error (aliquots measured to 4 decimal places below 30g, 2 decimal places for +30g)
 - Viscosity
 - Ability to adjust aliquots weighed in secondary container
- Reduce cost
 - Consumables – transfer pipettes vs pipette tips
 - Calibration of pipettes


Gravimetric Vs Volumetric Aliquot

Let's make a 100mL solution of 1000ppm Ca.
(Ca concentrate has concentration of 119,910ppm & density of 1.341 g/mL).

$$\left(\frac{100\text{mL} \times 1000\text{ppm}}{119910\text{ppm}} \right) = 0.83396\text{mL}$$

$$0.83396\text{mL} \times \frac{1.341\text{g}}{\text{mL}} = 1.1183\text{g}$$

Our volumetric aliquot is 0.8340mL – this should weigh 1.1183g



Theoretical aliquot =
1.1183g

From pipettor =
1.1107g

0.68% variation



Benefits of gravimetric preparation

- More reliable and consistent than volumetric preparation
 - Avoid pitfalls of plastic consumables
 - adhesion
 - Increased accuracy
 - Limit affect of atmospheric pressure and temperature
 - Small margin of error (aliquots measured to 4 decimal places below 30g, 2 decimal places for +30g)
 - Viscosity
 - Ability to adjust aliquots weighed in secondary container
- Reduce cost
 - Consumables – transfer pipettes vs pipette tips
 - Calibration of pipettes

Gravimetric Vs Volumetric Aliquot

Let's make a 100mL solution of 1000ppm Ca.
(Ca concentrate has concentration of 119,910ppm & density of 1.341 g/mL).

$$\left(\frac{100\text{mL} \times 1000\text{ppm}}{119910\text{ppm}} \right) = \mathbf{0.83396\text{mL}}$$

$$0.83396\text{mL} \times \frac{1.341\text{g}}{\text{mL}} = \mathbf{1.1183\text{g}}$$

Our volumetric aliquot is 0.8340mL – this should weigh 1.1183g



Theoretical aliquot =
1.1183g

From pipettor =
1.0744g

3.93% variation



Why buy custom multielement blends?



- Reliability
 - Traceable results – unbroken chain of comparisons
 - Documentation associated with CRMs – audits
 - Reduce interferences associated with stock MEBs
 - Typical uncertainty of 0.5%
 - ISO accreditations/registrations
- Save time and money
 - Reduce waste and number of solutions to purchase
 - Simplify working solution preparation
 - Avoid the need for in-house stability and self-assigned usage periods
 - **Transpiration Control Technology**




Scan me to submit a custom quote request!



inorganicventures.com

Transpiration Control Technology - TCT

- What is transpiration and why is it a concern?
- How does TCT solve that problem?
 - One-year vs five-year usage period
 - TCT allows more control over your method
 - Buy larger volume in series of smaller bottles
 - Maintain lot number
 - Reduce contamination potential
 - Control the expiration of your inventory



Scan me to navigate to our TCT education page!



[inorganicventures.com](https://www.inorganicventures.com)

Is TCT the only way to slow transpiration?

- Effects of transpiration can also be mitigated with refrigerated storage.
 - Allow standard to come to room temperature and mix by inversion 50x before use.
- Certificate of analysis also contains storage and handling information in Section 7.

- Store between approximately 4° - 30° C while in sealed TCT bag.

- While stored in the sealed TCT bag, transpiration of this CRM/RM is negligible. After opening the sealed TCT bag transpiration of the CRM/RM will occur, resulting in a gradual increase in the analyte concentration(s). It is the responsibility of the user to account for this effect. When the bottle is weighed both before and after being placed in storage, the mass difference observed will be a measure of transpiration mass loss.

- After opening the sealed TCT bag, keep cap tightly sealed when not in use and store between 4° - 24° C to minimize the effects of transpiration. Use at 20° ± 4° C to minimize volumetric dilution error when using the reported density. Do not pipette from the container. Do not return removed aliquots to container.



Visit our education
page for more
information and
resources on
frequently asked
questions!





Questions?



inorganicventures.com