

ICP Periodic Table Guide

1										2																																																																															
<div>H</div> <div>Hydrogen</div> <div>1.00794</div>										<div>He</div> <div>Helium</div> <div>4.002602</div>																																																																															
3		4												5		6		7		8		9		10																																																																	
<div>Li</div> <div>Lithium</div> <div>6.941</div>		<div>Be</div> <div>Beryllium</div> <div>9.012182</div>												<div>B</div> <div>Boron</div> <div>10.811</div>		<div>C</div> <div>Carbon</div> <div>12.0107</div>		<div>N</div> <div>Nitrogen</div> <div>14.0067</div>		<div>O</div> <div>Oxygen</div> <div>15.9994</div>		<div>F</div> <div>Fluorine</div> <div>18.9984032</div>		<div>Ne</div> <div>Neon</div> <div>20.1797</div>																																																																	
11										12										13										14										15										16										17										18																			
<div>Na</div> <div>Sodium</div> <div>22.989770</div>										<div>Mg</div> <div>Magnesium</div> <div>24.3050</div>																				<div>Al</div> <div>Aluminum</div> <div>26.981538</div>										<div>Si</div> <div>Silicon</div> <div>28.0855</div>										<div>P</div> <div>Phosphorus</div> <div>30.973761</div>										<div>S</div> <div>Sulfur</div> <div>32.065</div>										<div>Cl</div> <div>Chlorine</div> <div>35.453</div>										<div>Ar</div> <div>Argon</div> <div>39.948</div>									
19		20		21		22		23		24		25		26		27		28		29		30		31		32		33		34		35		36																																																							
<div>K</div> <div>Potassium</div> <div>39.0983</div>		<div>Ca</div> <div>Calcium</div> <div>40.078</div>		<div>Sc</div> <div>Scandium</div> <div>44.955910</div>		<div>Ti</div> <div>Titanium</div> <div>47.867</div>		<div>V</div> <div>Vanadium</div> <div>50.9415</div>		<div>Cr</div> <div>Chromium</div> <div>51.9961</div>		<div>Mn</div> <div>Manganese</div> <div>54.938049</div>		<div>Fe</div> <div>Iron</div> <div>55.845</div>		<div>Co</div> <div>Cobalt</div> <div>58.933200</div>		<div>Ni</div> <div>Nickel</div> <div>58.6934</div>		<div>Cu</div> <div>Copper</div> <div>63.546</div>		<div>Zn</div> <div>Zinc</div> <div>65.38</div>		<div>Ga</div> <div>Gallium</div> <div>69.723</div>		<div>Ge</div> <div>Germanium</div> <div>72.63</div>		<div>As</div> <div>Arsenic</div> <div>74.92160</div>		<div>Se</div> <div>Selenium</div> <div>78.971</div>		<div>Br</div> <div>Bromine</div> <div>79.904</div>		<div>Kr</div> <div>Krypton</div> <div>83.798</div>																																																							
37		38		39		40		41		42		43		44		45		46		47		48		49		50		51		52		53		54																																																							
<div>Rb</div> <div>Rubidium</div> <div>85.4678</div>		<div>Sr</div> <div>Strontium</div> <div>87.62</div>		<div>Y</div> <div>Yttrium</div> <div>88.90585</div>		<div>Zr</div> <div>Zirconium</div> <div>91.224</div>		<div>Nb</div> <div>Niobium</div> <div>92.90638</div>		<div>Mo</div> <div>Molybdenum</div> <div>95.96</div>		<div>Tc</div> <div>Technetium</div> <div>(98)</div>		<div>Ru</div> <div>Ruthenium</div> <div>101.07</div>		<div>Rh</div> <div>Rhodium</div> <div>102.90550</div>		<div>Pd</div> <div>Palladium</div> <div>106.42</div>		<div>Ag</div> <div>Silver</div> <div>107.8682</div>		<div>Cd</div> <div>Cadmium</div> <div>112.411</div>		<div>In</div> <div>Indium</div> <div>114.818</div>		<div>Sn</div> <div>Tin</div> <div>118.710</div>		<div>Sb</div> <div>Antimony</div> <div>121.760</div>		<div>Te</div> <div>Tellurium</div> <div>127.60</div>		<div>I</div> <div>Iodine</div> <div>126.90447</div>		<div>Xe</div> <div>Xenon</div> <div>131.293</div>																																																							
55		56		57-71		72		73		74		75		76		77		78		79		80		81		82		83		84		85		86																																																							
<div>Cs</div> <div>Cesium</div> <div>132.90545</div>		<div>Ba</div> <div>Barium</div> <div>137.327</div>		<div>La-Lu</div>		<div>Hf</div> <div>Hafnium</div> <div>178.49</div>		<div>Ta</div> <div>Tantalum</div> <div>180.9479</div>		<div>W</div> <div>Tungsten</div> <div>183.84</div>		<div>Re</div> <div>Rhenium</div> <div>186.207</div>		<div>Os</div> <div>Osmium</div> <div>190.23</div>		<div>Ir</div> <div>Iridium</div> <div>192.217</div>		<div>Pt</div> <div>Platinum</div> <div>195.078</div>		<div>Au</div> <div>Gold</div> <div>196.96655</div>		<div>Hg</div> <div>Mercury</div> <div>200.59</div>		<div>Tl</div> <div>Thallium</div> <div>204.3833</div>		<div>Pb</div> <div>Lead</div> <div>207.2</div>		<div>Bi</div> <div>Bismuth</div> <div>208.98038</div>		<div>Po</div> <div>Polonium</div> <div>(209)</div>		<div>At</div> <div>Astatine</div> <div>(210)</div>		<div>Rn</div> <div>Radon</div> <div>(222)</div>																																																							
87		88		89-103		104		105		106		107		108		109		110		111		112		113		114		115		116		117		118																																																							
<div>Fr</div> <div>Francium</div> <div>(223)</div>		<div>Ra</div> <div>Radium</div> <div>(226)</div>		<div>Ac-Lr</div>		<div>Rf</div> <div>Rutherfordium</div> <div>(261)</div>		<div>Db</div> <div>Dubnium</div> <div>(268)</div>		<div>Sg</div> <div>Seaborgium</div> <div>(269)</div>		<div>Bh</div> <div>Bohrium</div> <div>(270)</div>		<div>Hs</div> <div>Hassium</div> <div>(269)</div>		<div>Mt</div> <div>Mitnerium</div> <div>(276)</div>		<div>Ds</div> <div>Darmstadtium</div> <div>(281)</div>		<div>Rg</div> <div>Roentgenium</div> <div>(282)</div>		<div>Cn</div> <div>Coppernium</div> <div>(285)</div>		<div>Nh</div> <div>Nihonium</div> <div>(286)</div>		<div>Fl</div> <div>Flerovium</div> <div>(289)</div>		<div>Mc</div> <div>Moscovium</div> <div>(289)</div>		<div>Lv</div> <div>Livermorium</div> <div>(293)</div>		<div>Ts</div> <div>Tennessee</div> <div>(294)</div>		<div>Og</div> <div>Ognesson</div> <div>(294)</div>																																																							
57		58		59		60		61		62		63		64		65		66		67		68		69		70		71		72		73		74																																																							
<div>La</div> <div>Lanthanum</div> <div>138.9055</div>		<div>Ce</div> <div>Cerium</div> <div>140.116</div>		<div>Pr</div> <div>Praseodymium</div> <div>140.90765</div>		<div>Nd</div> <div>Neodymium</div> <div>144.24</div>		<div>Pm</div> <div>Promethium</div> <div>(145)</div>		<div>Sm</div> <div>Samarium</div> <div>150.36</div>		<div>Eu</div> <div>Europium</div> <div>151.964</div>		<div>Gd</div> <div>Gadolinium</div> <div>157.25</div>		<div>Tb</div> <div>Terbium</div> <div>158.92534</div>		<div>Dy</div> <div>Dysprosium</div> <div>162.500</div>		<div>Ho</div> <div>Holmium</div> <div>164.93032</div>		<div>Er</div> <div>Erbium</div> <div>167.259</div>		<div>Tm</div> <div>Thulium</div> <div>168.93421</div>		<div>Yb</div> <div>Ytterbium</div> <div>173.054</div>		<div>Lu</div> <div>Lutetium</div> <div>174.9668</div>																																																													
89		90		91		92		93		94		95		96		97		98		99		100		101		102		103		104		105																																																									
<div>Ac</div> <div>Actinium</div> <div>(227)</div>		<div>Th</div> <div>Thorium</div> <div>232.0381</div>		<div>Pa</div> <div>Protactinium</div> <div>231.03588</div>		<div>U</div> <div>Uranium</div> <div>238.02891</div>		<div>Np</div> <div>Neptunium</div> <div>(237)</div>		<div>Pu</div> <div>Plutonium</div> <div>(244)</div>		<div>Am</div> <div>Americium</div> <div>(243)</div>		<div>Cm</div> <div>Curium</div> <div>(247)</div>		<div>Bk</div> <div>Berkelium</div> <div>(247)</div>		<div>Cf</div> <div>Californium</div> <div>(251)</div>		<div>Es</div> <div>Einsteinium</div> <div>(252)</div>		<div>Fm</div> <div>Fermium</div> <div>(257)</div>		<div>Md</div> <div>Mendelevium</div> <div>(258)</div>		<div>No</div> <div>Nobelium</div> <div>(259)</div>		<div>Lr</div> <div>Lawrencium</div> <div>(260)</div>																																																													



Refine your results. Redefine your industry.

Inorganic Ventures

Since 1985, our team has perfected the process of manufacturing Certified Reference Materials, engineered to be stable, compatible, NIST traceable and are manufactured and tested under ISO 17034 & ISO 17025 guidelines (A2LA 883.01 & 883.02). We make almost any inorganic CRM for ICP, ICP-MS, IC, atomic absorption, wet chemistry, and QC applications, delivering ultimate standards across multiple industries.

We are chemists working for chemists to ensure the integrity of analytical testing, worldwide. Inorganic Ventures sets the standard for quality and accuracy. From our products and services to our educational commitment, you will not find a more dedicated analytical standards manufacturer in the world.

This guide includes essential data for 70+ elements, applicable to every ICP user. The analytical data presented includes chemical compatibilities, preferred emission lines, as well as major interferences and detection limits for both ICP-OES and ICP-MS. Learn more about solubility issues in different acid matrices, storage and handling tips, and the long-term stability of elements at different concentrations.

For a more interactive experience, visit our online Periodic Table at inorganicventures.com/periodic-table.



Contents

Li.....	4	Cd.....	40
Be.....	5	In.....	41
B.....	6	Sn.....	42
C.....	7	Sb.....	43
Na.....	8	Te.....	44
Mg.....	9	Cs.....	45
Al.....	10	Ba.....	46
Si.....	11	La.....	47
P.....	12	Hf.....	48
S.....	13	Ta.....	49
K.....	14	W.....	50
Ca.....	15	Re.....	51
Sc.....	16	Os.....	52
Ti.....	17	Ir.....	53
V.....	18	Pt.....	54
Cr.....	19	Au.....	55
Mn.....	20	Hg.....	56
Fe.....	21	Tl.....	57
Co.....	22	Pb.....	58
Ni.....	23	Bi.....	59
Cu.....	24	Ce.....	60
Zn.....	25	Pr.....	61
Ga.....	26	Nd.....	62
Ge.....	27	Sm.....	63
As.....	28	Eu.....	64
Se.....	29	Gd.....	65
Rb.....	30	Tb.....	66
Sr.....	31	Dy.....	67
Y.....	32	Ho.....	68
Zr.....	33	Er.....	69
Nb.....	34	Tm.....	70
Mo.....	35	Yb.....	71
Ru.....	36	Lu.....	72
Rh.....	37	Th.....	73
Pd.....	38	U.....	74
Ag.....	39		

Lithium

LOCATION: Group 1, Period 2

ATOMIC WEIGHT: 6.941

COORDINATION NUMBER: (6) (*coordination number in parentheses is assumed, not certain*)

CHEMICAL FORM IN SOLUTION: Li⁺(aq)
(large effective radius due to hydration sphere)

3	6.941
1342	1.0
180.7	
Li	
[He]2s	1
0.534	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO₃, H₂SO₄, and HF aqueous matrices. Stable with all metals and inorganic anions.

STABILITY: 2–100 ppb levels stable for months in 1% HNO₃/LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO₃/LDPE container.

Li CONTAINING SAMPLES (PREPARATION & SOLUTION):

Metal (dissolves very rapidly in water); Ores (sodium carbonate fusion in Pt⁰ followed by HCl dissolution – blank levels of Li in sodium carbonate critical); Organic Matrices (sulfuric/peroxide digestion or nitric/sulfuric/perchloric acid decomposition).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 670.784 nm	0.002/0.00002 µg/mL	1	atom	**
ICP-OES 460.286 nm	0.9/0.04 µg/mL	1	atom	Zr, Th
ICP-OES 323.261 nm	1.1/0.05 µg/mL	1	atom	<i>Sb, Th, Ni</i>
ICP-MS 7 amu	10 ppt	n/a	M+	

*ICP-OES D.L.'s are given as radial/axial view

**2nd order radiation from R.E.s on some optical designs

Beryllium

LOCATION: Group 2, Period 2

ATOMIC WEIGHT: 9.01218

COORDINATION NUMBER: 4

CHEMICAL FORM IN SOLUTION: $\text{Be}^+(\text{H}_2\text{O})_4^{+2}$

4

2472
1287

9.012
1.5

Be

[He]2s²
1.85

2

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , H_2SO_4 and HF aqueous matrices. Stable with all metals and inorganic anions.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 5–10% HNO_3 /LDPE container.

Be CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (is best dissolved in diluted H_2SO_4); BeO (boiling nitric, hydrochloric, or sulfuric acids or KHSO_4 fusion); Ores (H_2SO_4 /HF digestion or carbonate fusion in Pt^0); Organic Matrices (sulfuric/peroxide digestion or nitric/sulfuric/perchloric acid decomposition, or dry ash and dissolution according to the BeO procedure above).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 313.042 nm	0.0003/0.00009 $\mu\text{g/mL}$	1	ion	V, Ce, U
ICP-OES 234.861 nm	0.0003/0.00016 $\mu\text{g/mL}$	1	atom	Fe, Ta, Mo
ICP-OES 313.107 nm	0.0007/0.0005 $\mu\text{g/mL}$	1	ion	Ce, Th, Tm
ICP-MS 9 amu	4 ppt	n/a	M+	

*ICP-OES D.L.'s are given as radial/axial view

Boron

LOCATION: Group 13, Period 2

ATOMIC WEIGHT: 10.811

COORDINATION NUMBER: 4

CHEMICAL FORM IN SOLUTION: B(OH)_3 and B(OH)_4^{-1}

5	10.81
4002	2.0
2027	
B	
[He]2s ² 2p	
2.34	3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Moderately soluble in HCl , HNO_3 , H_2SO_4 and HF aqueous matrices and very soluble in NH_4OH . Stable with all metals and inorganic anions at low to moderate ppm levels.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–1,000 ppm solutions chemically stable for years in 1% HNO_3 /LDPE container. 1000–10,000 ppm stable for years in dilute NH_4OH /LDPE container.

B CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (crystalline form is scarcely attacked by acids or alkaline solutions; amorphous form is soluble in conc. HNO_3 or H_2SO_4); B(OH)_3 (water soluble); Ores (avoid acid digestions and use caustic fusions in Pt^0); Organic Matrices (dry ash mixed with Na_2CO_3 in Pt^0 at 450°C then increase heat to 1000°C to fuse; or perform a Na_2O_2 fusion in a Ni^0 crucible/Parr bomb).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 249.773 nm	0.003/0.001 $\mu\text{g/mL}$	1	atom	W, Ce, Co, Th, Ta, Mn, Mo, Fe
ICP-OES 249.678 nm	0.004/0.003 $\mu\text{g/mL}$	1	atom	Os, W, Co, Cr, Hf
ICP-OES 208.959 nm	0.007/0.0005 $\mu\text{g/mL}$	1	atom	Mo
ICP-MS 11 amu	700 ppt	n/a	M+	

*ICP-OES D.L.'s are given as radial/axial view

Carbon

LOCATION: Group 14, Period 2

ATOMIC WEIGHT: 12.011

COORDINATION NUMBER: 4

CHEMICAL FORM IN SOLUTION: (Carbon standard is made using Tartaric Acid)

6	12.011
4197	2.5
3827	
C	
[He]2s ² 2p ²	2,±4
2.25	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Moderately soluble in HCl, HNO₃, H₂SO₄ and HF aqueous matrices and very soluble in NH₄OH. Stable with all metals and inorganic anions at low to moderate ppm levels. Do not dilute or store Carbon standards using plastic containers or similar devices.

STABILITY: 2–100 ppb level stability unknown. 1000–10,000 ppm level stable for years in dilute acidic media in a glass container.

C CONTAINING SAMPLES (PREPARATION & SOLUTION): elemental amorphous or graphitic carbon (Oxidative closed vessel fusion such as a Na₂O₂ fusion in a sealed Ni⁰ crucible/Parr bomb); H₂CO₃ (water soluble); Organic Compounds (water solubility is best if possible, or perform a Na₂O₂ fusion in a sealed Ni⁰ crucible/Parr bomb).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 193.091 nm	0.05/0.005 µg/mL	1	atom	In, Ru, Mn
ICP-OES 247.856 nm	0.2/0.02 µg/mL	1	atom	Nb, V, Ti
ICP-MS 12 amu	(uncertain)	n/a	M+	

*ICP-OES D.L.'s are given as radial/axial view

Sodium

LOCATION: Group 1, Period 3

ATOMIC WEIGHT: 22.98977

COORDINATION NUMBER: (6) (*coordination number in parentheses is assumed, not certain*)

CHEMICAL FORM IN SOLUTION: Na⁺(aq) (largely ionic in nature)

11	22.990
883	1.0
98.0	
Na	
[Ne]3s	1
0.971	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO₃, H₂SO₄ and HF aqueous matrices. Stable with all metals and inorganic anions.

STABILITY: 2–100 ppb levels stable for months in 1% HNO₃/LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO₃/LDPE container.

Na CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (dissolves very rapidly in water); Ores (lithium carbonate fusion in graphite crucible followed by HCl dissolution – blank levels of Na in lithium carbonate critical); Organic Matrices (sulfuric/peroxide digestion or nitric/sulfuric/perchloric acid decomposition).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 589.595 nm	0.07/0.00009 µg/mL	1	atom	**
ICP-OES 588.995 nm	0.03/0.006 µg/mL	1	atom	**
ICP-OES 330.237 nm	2.0/0.09 µg/mL	1	atom	<i>Pd, Zn</i>
ICP-MS 23 amu	310 ppt	n/a	M+	⁴⁶ Ti ⁺² , ⁴⁶ Ca ⁺²

*ICP-OES D.L.'s are given as radial/axial view

**2nd order radiation from R.E.s on some optical designs

Magnesium

LOCATION: Group 2, Period 3

ATOMIC WEIGHT: 24.305

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Mg}(\text{H}_2\text{O})_6^{+2}$

12	24.305
1090	1.2
649	
Mg	
[Ne]2s ²	2
1.738	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , H_2SO_4 . Avoid HF, H_3PO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming insoluble silicates, carbonates, hydroxides, oxides, and tungstates in neutral and slightly acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–10% HNO_3 /LDPE container.

Mg CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (best dissolved in diluted HNO_3); Oxide (readily soluble in above compatible aqueous acidic solutions); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (sulfuric/peroxide digestion or nitric/sulfuric/perchloric acid decomposition, or dry ash and dissolution in dilute HCl. Do not heat when dissolving to avoid precipitation of SiO_2).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 279.553 nm	0.0002/0.00003 $\mu\text{g/mL}$	1	ion	Th
ICP-OES 280.270 nm	0.0003/0.00005 $\mu\text{g/mL}$	1	ion	U, V
ICP-OES 285.213 nm	0.002/0.00003 $\mu\text{g/mL}$	1	atom	U, Hf, Cr, Zr
ICP-MS 24 amu	42 ppt	n/a	M+	$^7\text{Li}^{17}\text{O}$, $^{48}\text{Ti}^{+2}$, $^{48}\text{Ca}^{+2}$

*ICP-OES D.L.'s are given as radial/axial view

Aluminum

LOCATION: Group 13, Period 3

ATOMIC WEIGHT: 26.98154

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Al}(\text{H}_2\text{O})_6^{+3}$

13	26.982
2520	1.5
660.25	
Al	
[Ne]3s23p ²	
2.699	3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , HF and H_2SO_4 . Avoid neutral media. Soluble in strongly basic NaOH forming the $\text{Al}(\text{OH})_4(\text{H}_2\text{O})_2^{1-}$ species. Stable with most metals and inorganic anions. The phosphate is insoluble in water and only slightly soluble in acid.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

AI CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (is best dissolved in HCl/ HNO_3); - Al_2O_3 (Na_2CO_3 fusion in Pt^0); - Al_2O_3 (soluble in acids such as HCl); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (sulfuric/peroxide digestion or nitric/sulfuric/perchloric acid decomposition, or dry ash and dissolution in dilute HCl).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 394.401 nm	0.05/0.006 $\mu\text{g/mL}$	1	atom	U, Ce
ICP-OES 396.152 nm	0.03/0.006 $\mu\text{g/mL}$	1	atom	Mo, Zr, Ce
ICP-OES 167.078 nm	0.1/0.009 $\mu\text{g/mL}$	1	ion	Fe
ICP-MS 27 amu	30 ppt	n/a	M+	¹² C ¹⁵ N, ¹³ C ¹⁴ N, ¹ H ¹² C ¹⁴ N, ¹¹ B ¹⁶ O, ⁵⁴ Cr ²⁺ , ⁵⁴ Fe ²⁺

*ICP-OES D.L.'s are given as radial/axial view

Silicon

LOCATION Group 14, Period 3

ATOMIC WEIGHT: 28.0855

COORDINATION NUMBER: 4

CHEMICAL FORM IN SOLUTION: $\text{Si}(\text{OH})_x(\text{F})_y^{2-}$

14

3267

1412

28.086

1.7

Si

[Ne]3s²3p²

2.33

4

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HF, H_3PO_4 , H_2SO_4 , and HNO_3 as the $\text{Si}(\text{OH})_x(\text{F})_y^{2-}$. Avoid neutral to basic media. Unstable at ppm levels with metals that would pull F- away (i.e. - do not mix with Alkaline or Rare Earths, or high levels of transition elements unless they are fluorinated). Stable with most inorganic anions with a tendency to hydrolyze forming silicic acid (silicic acid is soluble up to ~ 100 ppm in water) in all dilute acids except HF.

STABILITY: 2–100 ppb levels — stability unknown — (alone or mixed with all other metals) as the $\text{Si}(\text{OH})_x(\text{F})_y^{2-}$. 1–10,000 ppm single element solutions as the $\text{Si}(\text{OH})_x(\text{F})_y^{2-}$ chemically stable for years in 2–5% HNO_3 /trace HF in a LDPE container.

Si CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in 1:1:1 H_2O /HF/ HNO_3); Oxide — SiO_2 , amorphous (dissolve by heating in 1:1:1 H_2O /HF/ HNO_3); Oxide-quartz (fuse in Pt^0 with Na_2CO_3); Geological Samples (fuse in Pt^0 with Na_2CO_3 followed by HCl solution of the fuseate); Organic Matrices containing silicates and non volatile silicon compounds (dry ash at 450°C in Pt^0 and dissolve by gently warming with 1:1:1 H_2O /HF/ H_2SO_4 or fuse/ash with Na_2CO_3 and dissolve fuseate with HCl/ H_2O); Silicone Oils — dimethyl silicones depolymerize to form volatile monomer units when heated (measure directly in alcoholic KOH/xylene mixture where sample is treated first with the KOH at $60\text{--}100^\circ\text{C}$ to “unzip” the Si-O-Si polymeric structure or digest with conc. $\text{H}_2\text{SO}_4/\text{H}_2\text{O}_2$ followed by cooling and dissolution of the dehydrated silica with HF). Note that the direct analysis of silicone oils in an organic solvent will result in false high results due to high vapor pressure of volatile monomer units such as hexamethylcyclotrisiloxane. The KOH forms the $\text{K}_2+\text{Si}(\text{CH}_3)_2\text{O}=\text{NaCl}$, which is not volatile at room temperature.

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 251.611 nm	0.012/0.003 $\mu\text{g/mL}$	1	ion	Ta, U, Zn, Th
ICP-OES 212.412 nm	0.02/0.01 $\mu\text{g/mL}$	1	ion	Hf, Os, Mo, Ta
ICP-OES 288.158 nm	0.03/0.004 $\mu\text{g/mL}$	1	ion	Ta, Ce, Cr, Cd, Th
ICP-MS 28 amu	4000–8000 ppt	n/a	M+	$^{14}\text{N}_2$, $^{12}\text{C}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Phosphorus

LOCATION: Group 15, Period 3

ATOMIC WEIGHT: 30.97376

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: H_3PO_4

15
277
44.30
30.974
2.1

P

[Ne]3s²3p³
1.82
±3,4,5

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl , HNO_3 , H_2SO_4 , HF , water and NH_4OH . Stable with all metals and inorganic anions at low to moderate ppm levels under acidic conditions; precipitates with several metals occur in neutral media at higher concentrations.

STABILITY: 2–100 ppb levels — stability unknown — in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 0–1% HNO_3 /LDPE container.

P CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (never found free in nature); Oxides (water soluble); Ores (naturally occurring only as the phosphate, except for a few rare minerals found in meteorites — Na_2CO_3 fusion in Pt^0); Organic Matrices (dry ash mixed with Na_2CO_3 in Pt^0 at 450°C then increase heat to 1000°C to fuse; or, perform a $\text{H}_2\text{SO}_4/\text{H}_2\text{O}_2$ acid digestion).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 178.287 nm	0.03/0.002 $\mu\text{g/mL}$	1	atom	I
ICP-OES 177.495 nm	0.01/0.005 $\mu\text{g/mL}$	1	atom	Cu, Hf
ICP-OES 213.618 nm	0.08/0.03 $\mu\text{g/mL}$	1	atom	<i>Cu, Mo</i>
ICP-MS 31 amu	6000+ ppt	n/a	M+	¹⁵ N ₂ ¹ H, ¹⁵ N ¹⁶ O, ¹⁴ N ¹⁷ O, ¹³ C ¹⁸ O, ¹² C ¹⁸ O ¹ H, ⁶² Ni ²⁺

*ICP-OES D.L.'s are given as radial/axial view

Sulfur

LOCATION: Group 16, Period 3

ATOMIC WEIGHT: 32.066

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $(\text{O})_2\text{S}(\text{OH})_2$

16
32.07

444.75
2.4

115.36

S

[Ne]3s²3p⁴
±2,4,6

2.07

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO₃, H₃PO₄, and HF aqueous matrices, water, and NH₄OH. Stable with all metals and inorganic anions at low to moderate ppm levels under acidic conditions, except Ba, Pb, Ca, and to a lesser extent Sr.

STABILITY: 2–100 ppb levels — stability unknown — in 1% HNO₃/LDPE container. 1–10,000 ppm solutions chemically stable for years in LDPE container.

S CONTAINING SAMPLES (PREPARATION & SOLUTION): We most often get questions about the determination of S in rocks, silicates and insoluble sulfates (the finely powdered sample is fused in a Pt⁰ crucible with 6 times its weight of Na₂CO₃ + 0.5 grams KNO₃. The fuseate is extracted with water. Any BaSO₄ present in the sample is transposed by the carbonate fusion to the BaCO₃, which is left behind in the water-insoluble residue. If PbSO₄ is present, the fuseate should be boiled with a sodium carbonate saturated with CO₂ solution for 1 hour or more. The PbSO₄ will be transposed to the water insoluble carbonate which can be filtered off. Boiling the fuseate with a saturated carbonate solution is good insurance for samples containing Ba, Pb, Sr, and Ca. The Ba, Pb, Sr, and Ca free filtrate can be acidified and measured by ICP).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 166.669 nm	0.2/0.19 µg/mL	1	atom	Si, B
ICP-OES 182.034 nm	0.3/0.024 µg/mL	1	atom	
ICP-OES 143.328 nm	0.4/0.035 µg/mL	1	atom	
ICP-MS 32 amu	30,000 ppt	n/a	M+	¹⁶ O ₂ , ¹⁴ N ¹⁸ O, ¹⁵ N ¹⁷ O, ¹⁴ N ¹⁷ O ¹ H, ¹⁵ N ¹⁶ O ¹ H

*ICP-OES D.L.'s are given as radial/axial view

Potassium

LOCATION: Group 1, Period 4

ATOMIC WEIGHT: 39.0983

COORDINATION NUMBER: (6) (*coordination number in parentheses is assumed, not certain*)

CHEMICAL FORM IN SOLUTION: K+(aq)

19	39.098
759	0.9
63.35	
K	
[Ar]4s	1
0.86	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , H_2SO_4 , and HF aqueous matrices. Avoid use of HClO_4 due to insolubility of the perchlorate. Stable with all metals and inorganic anions except ClO_4^- .

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

K CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (dissolves very rapidly in water); Ores (sodium carbonate fusion in Pt^0 followed by HCl dissolution – blank levels of K in sodium carbonate critical); Organic Matrices (sulfuric/peroxide digestion).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 766.490 nm	0.4/0.001 $\mu\text{g/mL}$	1	atom	**
ICP-OES 771.531 nm	1.0/0.03 $\mu\text{g/mL}$	1	atom	**
ICP-OES 404.721 nm	1.1/0.05 $\mu\text{g/mL}$	1	atom	<i>U, Ce</i>
ICP-MS 39 amu	10 ppt	n/a	M+	$^{38}\text{Ar}^1\text{H}$, $^{23}\text{Na}^{16}\text{O}$, $^{78}\text{Se}^{+2}$
*ICP-OES D.L.'s are given as radial/axial view				

Calcium

LOCATION: Group 2, Period 4

ATOMIC WEIGHT: 40.078

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Ca}(\text{H}_2\text{O})_6^{+2}$

20	40.08
1484	1.0
839	
Ca	
[Ar]4s ²	2
1.55	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid H_2SO_4 , HF, H_3PO_4 and neutral to basic media. Stable with most metals and inorganic anions forming insoluble silicate, carbonate, hydroxide, oxide, fluoride, sulfate, oxalate, chromate, arsenate, and tungstate in neutral aqueous media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–10% HNO_3 /LDPE container.

Ca CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (best dissolved in diluted HNO_3); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolution in dilute HCl. Do not heat when dissolving to avoid precipitation of SiO_2). The oxide, hydroxide, carbonate, phosphate, and fluoride of calcium are soluble in % levels of HCl or HNO_3 . The sulfates (gypsum, anhydrite, etc.), certain silicates, and complex compounds require fusion with Na_2CO_3 followed by HCl/water dissolution. Note that contamination is a very real problem when analyzing for trace levels.

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 393.366 nm	0.0002/0.00004 $\mu\text{g/mL}$	1	ion	U, Ce
ICP-OES 396.847 nm	0.0005/0.00006 $\mu\text{g/mL}$	1	ion	Th
ICP-OES 422.673 nm	0.01/0.001 $\mu\text{g/mL}$	1	atom	Ge
ICP-MS 44 amu	1200 ppt	n/a	M+	$^{16}\text{O}_2^{12}\text{C}$, $^{28}\text{Si}^{16}\text{O}$, $^{88}\text{Sr}^{+2}$

*ICP-OES D.L.'s are given as radial/axial view

Scandium

LOCATION: Group 3, Period 4

ATOMIC WEIGHT: 44.95591

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Sc}(\text{H}_2\text{O})_6^{+2}$

21
2831
1539
44.956
1.3

Sc

[Ar]3d4s²
2.99
3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl , H_2SO_4 , and HNO_3 . Avoid HF , H_3PO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride. The fluoride is soluble in excess HF , forming ScF_6^{3-} (not recommended for standard preparations).

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 5–10% HNO_3 /LDPE container. Small atomic radius increases hydrolysis requiring higher acid levels than other Rare Earths.

Sc CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3 – aqua regia or nitric/perchloric/sulfuric acid digestions can be used – exercise caution when using perchloric acid).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 335.373 nm	0.004/0.00002 $\mu\text{g/mL}$	1	ion	
ICP-OES 337.215 nm	0.004/0.00002 $\mu\text{g/mL}$	1	ion	Ti, U, Ni, Rh
ICP-OES 424.683 nm	0.003/0.00002 $\mu\text{g/mL}$	1	ion	Ce
ICP-MS 45 amu	2.3 ppt	n/a	M+	$^{16}\text{O}_2^{12}\text{C}^{1}\text{H}$, $^{29}\text{Si}^{16}\text{O}$, $^{90}\text{Zr}^{+2}$

*ICP-OES D.L.'s are given as radial/axial view

Titanium

LOCATION: Group 4, Period 4

ATOMIC WEIGHT: 47.867

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Ti}(\text{F})_6^{-2}$

22
3289
1670

47.867

Ti

[Ar]3d²4s²
4.50

3,4

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in concentrated HCl, HF, H₃PO₄, H₂SO₄, and HNO₃. Avoid neutral to basic media. Unstable at ppm levels with metals that would pull F- away (i.e. do not mix with Alkaline or Rare Earths or high levels of transition elements unless they are fluorinated). Stable with most inorganic anions with a tendency to hydrolyze forming the hydrated oxide in all dilute acids except HF.

STABILITY: 2–100 ppb levels stable (alone or mixed with all other metals) as the $\text{Ti}(\text{F})_6^{-2}$ for months in 1% HNO₃/LDPE container. 1–10,000 ppm single element solutions as the $\text{Ti}(\text{F})_6^{-2}$ chemically stable for years in 2–5% HNO₃/trace HF in an LDPE container.

Ti CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in H₂O/HF CAUTION — *powder reacts violently*); Oxide — low temperature history *anatase* or *rutile* (dissolved by heating in 1:1:1 H₂O/HF/H₂SO₄); Oxide — high temperature history {~ 800°C} *brookite* (fuse in Pt⁰ with K₂S₂O₇); Ores (fuse in Pt⁰ with KF + K₂S₂O₇ — no KF if silica not present); Organic Matrices (dry ash at 450°C in Pt⁰ and dissolve by heating with 1:1:1 H₂O/HF/H₂SO₄ or fuse ash with pyrosulfate if oxide is as plastic pigment and likely in brookite crystalline form).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 334.941 nm	0.0038/0.000028 µg/mL	1	ion	Nb, Ta, Cr, U
ICP-OES 336.121 nm	0.0053/0.000034 µg/mL	1	ion	W, Mo, Co
ICP-OES 323.452 nm	0.0054/0.00092 µg/mL	1	ion	Ce, Ar, Ni
ICP-MS 48 amu	14 ppt	n/a	M+	³² S ¹⁶ O, ³⁴ S ¹⁴ N, ¹⁴ N ¹⁶ O ¹⁸ O, ¹⁴ N ¹⁷ O ₂ , ³⁶ Ar ¹² C, ⁴⁸ Ca, [⁹⁶ X=2 (where X = Zr, Mo, Ru)]

*ICP-OES D.L.'s are given as radial/axial view

Vanadium

LOCATION: Group 5, Period 4

ATOMIC WEIGHT: 50.9416

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{H}_2\text{V}_{10}\text{O}_{28}^{4-}$

23
3409
1902

50.942
1.5

V

[Ar]3d³4s²
6.11

2,3,4,5

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , H_2SO_4 , HF, H_3PO_4 , and strong basic media. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

V CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (fusion with NaOH or KOH in Ni^0 or $\text{Na}_2\text{CO}_3/\text{KNO}_3$); Oxides (V_2O_3 – use HCl; V_2O_4 – use HCl or HNO_3 ; V_2O_5 - use conc. acids); Ores (Na_2CO_3 / KNO_3 in Pt^0 (*caution – nitrates attack Pt^0*) followed by water extraction of fuseate); Organic Matrices (ash at 450°C followed by dissolving according to V_2O_5 above).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 309.311 nm	0.005/0.001 $\mu\text{g/mL}$	1	ion	Mg, U, Th
ICP-OES 292.402 nm	0.006/0.001 $\mu\text{g/mL}$	1	ion	Th
ICP-OES 290.882 nm	0.008/0.0008 $\mu\text{g/mL}$	1	atom	Hf, Nb
ICP-MS 51 amu	4 ppt	n/a	M+	$^{34}\text{S}^{16}\text{O}^{16}\text{H}$, $^{35}\text{Cl}^{16}\text{O}$, $^{38}\text{Ar}^{13}\text{C}$, $^{36}\text{Ar}^{15}\text{N}$, $^{36}\text{Ar}^{14}\text{N}^{16}\text{H}$, $^{37}\text{Cl}^{14}\text{N}$, $^{36}\text{S}^{15}\text{N}$, $^{33}\text{S}^{18}\text{O}$, $^{34}\text{S}^{17}\text{O}$, $^{102}\text{Ru}^{+2}$, $^{102}\text{Pd}^{+2}$

*ICP-OES D.L.'s are given as radial/axial view

Chromium

LOCATION: Group 6, Period 4

ATOMIC WEIGHT: 51.9961

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Cr}(\text{H}_2\text{O})_6^{3+}$

24	51.996
2672	1.6
1857	
Cr	
[Ar]3d ⁵ 4s	2,3,6
7.19	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl , HNO_3 , H_2SO_4 , HF , H_3PO_4 . Avoid basic media. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1-5% HNO_3 /LDPE container.

Cr CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HCl); Oxides / Ores (chrome ore/oxides are very difficult to dissolve). The following procedures A–D are commonly used:

- Fusion with KHSO_4 and extraction with hot KCl . The residue fused with Na_2CO_3 and KClO_3 , 3:1
- Fusion with NaKSO_4 and NaF , 2:1
- Fusion with magnesite or lime and sodium or potassium carbonates, 4:1
- Fusion with Na_2O_2 or NaOH and KNO_3 or NaOH and Na_2O_2 .

Nickel, iron, copper, or silver crucibles should be used for D. Platinum may be used for A, B, and C; Organic Matrices (ash at 450°C followed by one of the fusion methods above or sulfuric/hydrogen peroxide acid digestions *may* be applicable to non-oxide containing samples).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 205.552 nm	0.006/0.0008 $\mu\text{g/mL}$	1	ion	Os
ICP-OES 284.325 nm	0.008/0.0007 $\mu\text{g/mL}$	1	ion	
ICP-OES 276.654 nm	0.01/0.001 $\mu\text{g/mL}$	1	ion	Cu, Ta, V
ICP-MS 52 amu	40 ppt	n/a	M-	$^{36}\text{S}^{16}\text{O}$, $^{36}\text{Ar}^{16}\text{O}^{**}$

*ICP-OES D.L.'s are given as radial/axial view

**The 50Cr, 53Cr, 54Cr lines suffer from many more potential interferences from sulfur, chlorine, and argon compounds of oxygen, nitrogen, and carbon.

Manganese

LOCATION: Group 7, Period 4

ATOMIC WEIGHT: 54.9380

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Mn}(\text{H}_2\text{O})_6^{2+}$

25
2062
1244
54.938
1.6

Mn

[Ar]3d⁵4s²
7.43
2,3,4,6,7

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , H_2SO_4 , HF, H_3PO_4 . Avoid basic media. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Mn CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in dilute acids); Oxides (soluble in dilute acids); Ores (dissolve with HCl. If silica is present, add HF and then fume off silica by adding H_2SO_4 and heat to SO_3 fumes – dense white fumes).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 257.610 nm	0.0014/0.00002 µg/mL	1	ion	Ce, W, Re
ICP-OES 259.373 nm	0.0016/0.00002 µg/mL	1	ion	U, Ta, Mo, Fe, Nb
ICP-OES 260.569 nm	0.0021/0.00002 µg/mL	1	ion	Co
ICP-MS 55 amu	10 ppt	n/a	M+	⁴⁰ Ar ¹⁴ N ¹ H, ³⁹ K ¹⁶ O, ³⁷ Cl ¹⁸ O, ⁴⁰ Ar ¹⁵ N, ³⁸ Ar ¹⁷ O, ³⁶ Ar ¹⁸ O ¹ H, ³⁸ Ar ¹⁶ O ¹ H, ³⁷ Cl ¹⁷ O ¹ H, ²³ Na ³² S
*ICP-OES D.L.'s are given as radial/axial view				

Iron

LOCATION: Group 8, Period 4

ATOMIC WEIGHT: 55.847

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Fe}(\text{H}_2\text{O})_6^{3+}$

26	55.847
2862	1.6
1563	
Fe	
[Ar]3d ⁶ 4s ²	2,3
7.86	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl , HNO_3 , H_2SO_4 , HF , H_3PO_4 . Avoid basic media. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Fe CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HCl); Oxides (if the oxide has been at a high temperature then Na_2CO_3 fusion in Pt^0 followed by HCl dissolution, otherwise dissolve in dilute HCl); Ores (see Oxides above using only the fusion approach).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 238.204 nm	0.005/0.001 $\mu\text{g/mL}$	1	ion	Ru, Co
ICP-OES 239.562 nm	0.005/0.001 $\mu\text{g/mL}$	1	ion	Co, W, Cr
ICP-OES 259.940 nm	0.006/0.001 $\mu\text{g/mL}$	1	ion	Hf, Nb
ICP-MS 56 amu	970 ppt	n/a	M+	⁴⁰ Ar ¹⁵ N ¹ H, ⁴⁰ Ar ¹⁶ O, ³⁸ Ar ¹⁸ O, ³⁷ Cl ¹⁸ O ¹ H, ⁴⁰ Ca ¹⁶ O

*ICP-OES D.L.'s are given as radial/axial view

Cobalt

LOCATION: Group 9, Period 4

ATOMIC WEIGHT: 58.9332

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Co}(\text{H}_2\text{O})_6^{2+}$

27
2928
1495
58.933
1.7

Co

[Ar]3d⁷4s²
8.9
2,3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , H_2SO_4 , HF, H_3PO_4 . Avoid basic media. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Co CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); Oxides (soluble in HCl); Ores (dissolve in HCl/ HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 238.892 nm	0.01/0.002 µg/mL	1	ion	Fe, W, Ta
ICP-OES 228.616 nm	0.01/0.001 µg/mL	1	ion	
ICP-OES 237.862 nm	0.01/0.002 µg/mL	1	ion	W, Re, Al, Ta
ICP-MS 59 amu	2 ppt	n/a	M+	⁴² Ca ¹⁶ O ¹ H, ⁴⁰ Ar ¹⁸ O ¹ H, ³⁶ Ar ²³ Na, ⁴³ Ca ¹⁶ O, ²⁴ Mg ³⁵ Cl
*ICP-OES D.L.'s are given as radial/axial view				

Nickel

LOCATION: Group 10, Period 4

ATOMIC WEIGHT: 58.6934

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Ni}(\text{H}_2\text{O})_6^{2+}$

28
2914
1453
58.6934
1.8

Ni

[Ar]3d⁸4s²
8.9
2,3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl , HNO_3 , H_2SO_4 , HF , H_3PO_4 . Avoid basic media. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Ni CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); Oxides (soluble in HCl); Ores (dissolve in HCl/HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 221.647 nm	0.01/0.0009 $\mu\text{g}/\text{mL}$	1	ion	Si
ICP-OES 232.003 nm	0.02/0.006 $\mu\text{g}/\text{mL}$	1	atom	Cr, Re, Os, Nb, Ag, Pt, Fe
ICP-OES 231.604 nm	0.02/0.002 $\mu\text{g}/\text{mL}$	1	ion	Sb, Ta, Co
ICP-MS 60 amu	100 ppt	n/a	M+	⁴³ Ca ¹⁶ O ¹ H, ⁴⁴ Ca ¹⁶ O, ²³ Na ³⁷ Cl

*ICP-OES D.L.'s are given as radial/axial view

Copper

LOCATION: Group 11, Period 4

ATOMIC WEIGHT: 63.546

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Cu}(\text{H}_2\text{O})_6^{2+}$

29
2563
1084.6
63.546
1.8

Cu

[Ar]3d¹⁰4s
8.96
1,2

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl , HNO_3 , H_2SO_4 , HF , H_3PO_4 . Avoid basic media. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Cu CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); Oxides (soluble in HCl); Ores (dissolve in HCl/HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text* indicates severe at ~ concs.)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 324.754 nm	0.06/0.001 $\mu\text{g}/\text{mL}$	1	atom	Nb, U, Th, Mo, Hf
ICP-OES 224.700 nm	0.01/0.001 $\mu\text{g}/\text{mL}$	1	ion	<i>Pb</i> , Ir, Ni, W
ICP-OES 219.958 nm	0.01/0.002 $\mu\text{g}/\text{mL}$	1	atom	Th, Ta, Nb, U, Hf
ICP-MS 63 amu	10 ppt	n/a	M+	⁴⁰ Ar ²³ Na, ⁴⁷ Ti ¹⁶ O, ¹⁴ N ¹² C ³⁷ Cl, ¹⁶ O ¹² C- ³⁵ Cl, ⁴⁴ Ca ¹⁸ O ¹ H, ²³ Na ⁴⁰ Ca
*ICP-OES D.L.'s are given as radial/axial view				

Zinc

LOCATION: Group 12, Period 4

ATOMIC WEIGHT: 65.38

COORDINATION NUMBER: 4

CHEMICAL FORM IN SOLUTION: Zn(OH)(aq)^{1+}

30
907
419.73
65.38
1.7

Zn

[Ar]3d¹⁰4s²
7.13
2

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl , HNO_3 , H_2SO_4 , HF , H_3PO_4 . Avoid basic media that promotes formation of insoluble carbonate and hydroxide. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1-5% HNO_3 /LDPE container.

Zn CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); Oxides (soluble in HCl); Ores (dissolve in HCl/HNO_3); Organic based (dry ash at 450°C and dissolve ash in HCl), (sulfuric/peroxide acid digestion).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 213.856 nm	0.002/0.0004 $\mu\text{g/mL}$	1	atom	Ni, Cu, V
ICP-OES 202.548 nm	0.004/0.0002 $\mu\text{g/mL}$	1	ion	Nb, Cu, Co, Hf
ICP-OES 206.200 nm	0.006/0.0006 $\mu\text{g/mL}$	1	ion	Sb, Ta, Bi, Os
ICP-MS 66 amu	7 ppt	n/a	M-	⁵⁰ Ti ¹⁶ O, ⁵⁰ Cr ¹⁶ O, ⁵⁰ V ¹⁶ O, ³⁴ S ¹⁶ O ₂ , ³² S ¹⁶ O ¹⁸ O, ³² S ¹⁷ O ₂ , ³³ S ¹⁶ O ¹⁷ O, ³² S ³⁴ S, ³³ S ₂

*ICP-OES D.L.'s are given as radial/axial view

Gallium

LOCATION: Group 13, Period 4

ATOMIC WEIGHT: 69.723

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Ga}(\text{H}_2\text{O})_6^{+3}$

31	69.72
2205	1.8
29.9	
Ga	
[Ar]3d ¹⁰ 4s ² 4p	
5.904	2

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , and H_2SO_4 . Avoid neutral media. Stable with most metals and inorganic anions. The fluoride is insoluble in water but soluble in HF.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Ga CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (is best dissolved in HCl/ HNO_3); Ga_2O_3 (Na_2CO_3 fusion in Pt^0); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (sulfuric/peroxide digestion or nitric/sulfuric/perchloric acid decomposition or dry ash and dissolution in dilute HCl).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 141.444 nm	0.05/0.001 $\mu\text{g}/\text{mL}$	1	ion	Hg
ICP-OES 294.364 nm	0.05/0.008 $\mu\text{g}/\text{mL}$	1	atom	Ce, U, Ni
ICP-OES 417.206 nm	0.07/0.005 $\mu\text{g}/\text{mL}$	1	atom	Ti, Ce
ICP-MS 69 amu	2 ppt	n/a	M+	$^{35}\text{Cl}^{16}\text{O}^{18}\text{O}$, $^{35}\text{Cl}^{17}\text{O}_2$, $^{37}\text{Cl}^{16}\text{O}_2$, $^{36}\text{Ar}^{33}\text{S}$, $^{33}\text{S}^{18}\text{O}_2$, $^{34}\text{S}^{17}\text{O}^{18}\text{O}$, $^{36}\text{S}^{16}\text{O}^{17}\text{O}$, $^{33}\text{S}^{36}\text{S}$, $^{53}\text{Cr}^{16}\text{O}$, [$^{138}\text{X}^{2+}$ (where X = Ba, La, Ce)]

*ICP-OES D.L.'s are given as radial/axial view

Germanium

LOCATION: Group 14, Period 4

ATOMIC WEIGHT: 72.630

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Ge}(\text{OH})_x(\text{F})_y^{2-}$

32	72.63
2834	2.0
937.4	
Ge	
[Ar]3d ¹⁰ 4s ² 4p ²	
5.32	4

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl, HF, H_3PO_4 , H_2SO_4 , and HNO_3 as the $\text{Ge}(\text{OH})_x(\text{F})_y^{2-}$. Avoid neutral to basic media. Unstable at ppm levels with metals that would pull F- away (i.e. Do not mix with Alkaline or Rare Earths or high levels of transition elements unless they are fluorinated). Stable with most inorganic anions with a tendency to hydrolyze.

STABILITY: 2–100 ppb levels — stability unknown alone or mixed with all other metals as the $\text{Ge}(\text{OH})_x(\text{F})_y^{2-}$. 1–10,000 ppm single element solutions as the $\text{Ge}(\text{OH})_x(\text{F})_y^{2-}$ chemically stable for years in 2–5% HNO_3 / trace HF in a LDPE container.

Ge CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in 1:1:1 $\text{H}_2\text{O}/\text{HF}/\text{HNO}_3$); Oxide — GeO (readily soluble in HCl or NaOH), GeO_2 (fuse in Pt^0 with Na_2CO_3 followed by HCl solution of the fuseate); Geological Samples (fuse in Pt^0 with Na_2CO_3 followed by HCl solution of the fuseate); Organic Matrices (dry ash at 450°C in Pt^0 and dissolve by gently warming with 1:1:1 $\text{H}_2\text{O}/\text{HF}/\text{H}_2\text{SO}_4$ or fuse ash with Na_2CO_3 and dissolve fuseate with HCl/ H_2O).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 164.919 nm	0.01/0.001 µg/mL	1	ion	Co, Fe, Cu
ICP-OES 219.871 nm	0.06/0.009 µg/mL	1	atom	W, Ir, Re, Co
ICP-OES 265.117 nm	0.05/0.009 µg/mL	1	atom	Ir, Re
ICP-MS 72 amu	20 ppt	n/a	M+	³⁶ Ar ₂ , ³⁷ Cl ¹⁷ O ¹⁸ O, ³⁷ Cl ³⁵ Cl, ³⁶ S ¹⁸ O ₂ , ³⁶ S ₂ , ³⁶ Ar ³⁶ S, ⁵⁶ Fe ¹⁶ O, ⁴⁰ Ar ¹⁶ O ₂ , ⁴⁰ Ca ¹⁶ O ₂ , ⁴⁰ Ar ³² S, ¹⁴⁴ Nd ²⁺ , ¹⁴⁴ Sm ²⁺

*ICP-OES D.L.'s are given as radial/axial view

Arsenic

LOCATION: Group 15, Period 4

ATOMIC WEIGHT: 74.9216

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: H_3AsO_4 and HAsO_2

33
603 (subl.)
808 (28 atm)

74.922
2.2

As

[Ar]3d¹⁰4s²4p³
5.73

±3,5

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Arsenic has no cationic chemistry. It is soluble in HCl , HNO_3 , H_3PO_4 , H_2SO_4 and HF aqueous matrices water and NH_4OH . It is stable with most inorganic anions (forms arsenate when boiled with chromate) but many cationic metals form the insoluble arsenates under pH neutral conditions. When fluorinated and/or under acidic conditions arsenate formation is typically not a problem at moderate to low concentrations.

STABILITY: 2–100 ppb levels — stable for months alone or mixed with other elements at equivalent levels — in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1-5% HNO_3 /LDPE container.

As CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in 1:1 $\text{H}_2\text{O}/\text{HNO}_3$); Oxides (the oxide exists in crystalline and amorphous forms where the amorphous form is more water soluble. The oxides typically dissolve in dilute acidic solutions when boiled); Minerals (one gram of powdered sample is fused in a Ni^0 crucible with 10 grams of a 1:1 mix of K_2CO_3 and KNO_3 and the melt extracted with hot water); Organic Matrices (0.2 to 0.5 grams of the sample are fused with 15 grams of a 1:1 $\text{Na}_2\text{CO}_3/\text{Na}_2\text{O}_2$ mix in a Ni^0 crucible. The fuseate is extracted with water and acidified with HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text* indicates severe at ~ concs.)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 189.042 nm	0.05/0.005 $\mu\text{g}/\text{mL}$	1	atom	Cr
ICP-OES 193.696 nm	0.1/0.01 $\mu\text{g}/\text{mL}$	1	atom	V, Ge
IICP-OES 228.812 nm	0.1/0.01 $\mu\text{g}/\text{mL}$	1	atom	Cd, Pt, Ir, Co
IICP-MS 75 amu	30 ppt	n/a	M+	⁴⁰ Ar ³⁵ Cl, ⁵⁹ Co ¹⁶ O, ³⁶ Ar ³⁸ Ar ¹ H, ³⁸ Ar ³⁷ Cl, ³⁶ Ar ³⁹ K, ¹⁵⁰ Nd ²⁺ , ¹⁵⁰ Sm ²⁺

*ICP-OES D.L.'s are given as radial/axial view

Selenium

LOCATION: Group 16, Period 4

ATOMIC WEIGHT: 78.971

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: H_2SeO_3

34	78.97
685	2.5
221	
Se	
[Ar]3d ¹⁰ 4s ² 4p ² 4.79 -2,4,6	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl , HNO_3 , H_3PO_4 , H_2SO_4 and HF aqueous matrices and water. It is stable with most inorganic anions but many cationic metals form the insoluble selenites under pH neutral conditions. When fluorinated and/or under acidic conditions precipitation is typically not a problem at moderate to low concentrations.

STABILITY: 2–100 ppb levels stable for months alone or mixed with other elements at equivalent levels in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Se CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); Oxides (readily soluble in water); Minerals and alloys (acid digestion with HNO_3 or HNO_3/HF); Organic Matrices (acid digestion with hot concentrated H_2SO_4 accompanied by the careful dropwise addition of H_2O_2 until clear).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 196.026 nm	0.08/0.006 $\mu\text{g/mL}$	1	atom	Fe
ICP-OES 203.985 nm	0.2/0.05 $\mu\text{g/mL}$	1	atom	<i>Sb, Ir, Cr, Ta</i>
ICP-OES 206.279 nm	0.3/0.16 $\mu\text{g/mL}$	1	atom	<i>Cr, Pt</i>
ICP-MS 82 amu	200 ppt	n/a	M+	¹² C ³⁵ Cl ₂

*ICP-OES D.L.'s are given as radial/axial view

Rubidium

LOCATION: Group 1, Period 5

ATOMIC WEIGHT: 85.4678

COORDINATION NUMBER: (6) (*coordination number in parentheses is assumed, not certain*)

CHEMICAL FORM IN SOLUTION: Rb(aq)

37	85.468
688	0.9
39.64	
Rb	
[Kr]5s	
1.532	1

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , H_2SO_4 , and HF aqueous matrices. Stable with most metals and inorganic anions. Forms insoluble $\text{Rb}_2[\text{PtCl}_6]$ (0.028g/100mL ^{20}aq).

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 0.1–1% HNO_3 /LDPE container.

Rb CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (dissolves very rapidly in water); Ores (sodium carbonate fusion in Pt^0 followed by HCl dissolution – blank levels of Rb in sodium carbonate critical); Organic Matrices (sulfuric/peroxide digestion or nitric/sulfuric/perchloric acid decomposition).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 420.185 nm	40/10 $\mu\text{g/mL}$	1	atom	Fe, Zr
ICP-MS 85 amu	1.5 ppt	n/a	M+	$^{69}\text{Ga}^{16}\text{O}$, $^{170}\text{Er}^{+2}$, $^{170}\text{Yb}^{+2}$
*ICP-OES D.L.'s are given as radial/axial view				

Strontium

LOCATION: Group 2, Period 5

ATOMIC WEIGHT: 87.62

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Sr}(\text{H}_2\text{O})_6^{+2}$

38	87.62
1377	1.0
768	
Sr	
[Kr]5s ²	2
2.54	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid H_2SO_4 , HF and neutral to basic media. Stable with most metals and inorganic anions forming insoluble silicate, carbonate, hydroxide, oxide, fluoride, sulfate, oxalate, chromate, arsenate, and tungstate in neutral aqueous media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Sr CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (is best dissolved in diluted HNO_3); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolution in dilute HCl. Do not heat when dissolving to avoid precipitation of SiO_2).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 407.771 nm	0.0004/0.00006 $\mu\text{g/mL}$	1	ion	U, Ce
ICP-OES 421.552 nm	0.0008/0.00004 $\mu\text{g/mL}$	1	ion	Rb
ICP-OES 460.733 nm	0.07/0.003 $\mu\text{g/mL}$	1	atom	Ce
ICP-MS 88 amu	1200 ppt	n/a	M+	⁷² Ge ¹⁶ O, ¹⁷⁶ Yb ⁺² , ¹⁷⁶ Lu ⁺² , ¹⁷⁶ Hf ⁺²

*ICP-OES D.L.'s are given as radial/axial view

Yttrium

LOCATION: Group 3, Period 5

ATOMIC WEIGHT: 88.906

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Y(OH)(H}_2\text{O)}_x^{+2}$

39	88.906
3338	1.1
1526	
Y	
[Kr]4d5s ²	
4.47	3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl , H_2SO_4 , and HNO_3 . Avoid HF , H_3PO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Y CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 360.073 nm	0.005/0.000036 $\mu\text{g/mL}$	1	ion	Ce, Th
ICP-OES 371.030 nm	0.004/0.00007 $\mu\text{g/mL}$	1	ion	Ce
ICP-OES 377.433 nm	0.005/0.0009 $\mu\text{g/mL}$	1	ion	Ta, Th
ICP-MS 89 amu	0.8 ppt	n/a	M+	$^{73}\text{Ge}^{16}\text{O}$, $^{178}\text{Hf}^{+2}$

*ICP-OES D.L.'s are given as radial/axial view

Zirconium

LOCATION: Group 4, Period 5

ATOMIC WEIGHT: 91.224

COORDINATION NUMBER: 6, 7, 8 (*coordination numbers 7, 8 are observed less frequently*)

CHEMICAL FORM IN SOLUTION: $\text{Zr}(\text{F})_6^{-2}$

40	91.22
4409	1.2
1852	
Zr	
[Kr]4d ² 5s ²	
6.51	4

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in concentrated HCl, HF, H_2SO_4 (very hot) and HNO_3 . Avoid H_3PO_4 and neutral to basic media. Unstable at ppm levels with metals that would pull F⁻ away (i.e. do not mix with Alkaline or Rare Earths or high levels of transition elements unless they are fluorinated). Stable with most inorganic anions but precipitation with phosphate, oxalate, and tartrate with a tendency to hydrolyze forming the hydrated oxide in all dilute acids except HF.

STABILITY: 2–100 ppb levels stable (alone or mixed with all other metals that are at comparable levels) as the $\text{Zr}(\text{F})_6^{-2} + \text{Zr}(\text{OH})_4\text{F}_2^{-2}$ for months in 1% HNO_3 /LDPE container. 1–10,000 ppm single element solutions as the $\text{Zr}(\text{F})_6^{-2}$ chemically stable for years in 2–5% HNO_3 /trace HF in an LDPE container.

Zr CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in H_2O /HF/ HNO_3); Oxide unlike TiO_2 , the ZrO_2 is best fused in one of the following ways (Na_2O_2 in Ni^0 , Na_2CO_3 in Pt^0 or Borax in Pt^0); Organic Matrices (dry ash at 450°C in Pt^0 and dissolve by fusing with Na_2CO_3 and dissolving in HF/ HNO_3 / H_2O).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 343.823 nm	0.007/0.0004 $\mu\text{g/mL}$	1	ion	Hf, Nb
ICP-OES 339.198 nm	0.008/0.0007 $\mu\text{g/mL}$	1	ion	<i>Th</i> , Mo
ICP-OES 272.261 nm	0.018/0.001 $\mu\text{g/mL}$	1	ion	<i>Cr</i> , V, <i>Th</i> , W
ICP-MS 90 amu	2 ppt	n/a	M+	$^{74}\text{Ge}^{16}\text{O}$, $^{74}\text{Se}^{16}\text{O}$, [$^{180}\text{X}^{+2}$ (where X = Hf, Ta, W)]

*ICP-OES D.L.'s are given as radial/axial view

Niobium

LOCATION: Group 5, Period 5

ATOMIC WEIGHT: 92.9064

COORDINATION NUMBER: 6, 7, 8 (*coordination numbers 7, 8 are observed less frequently*)

CHEMICAL FORM IN SOLUTION: NbOF_5^{-2}

41	92.906
4744	1.2
2467	
Nb	
[Kr]4d ⁴ 5s	3,5
8.57	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in concentrated HCl and dilute HF/HNO₃. Avoid neutral to basic media. Unstable at ppm levels with metals that would pull F- away (i.e. do not mix with Alkaline or Rare Earths or high levels of transition elements unless they are fluorinated). Stable with most inorganic anions provided it is in the chemical form shown above.

STABILITY: 2–100 ppb levels stable (alone or mixed with all other metals that are at comparable levels) as the NbOF_5^{-2} for 5 months in 1% HNO₃/LDPE container. 1–10,000 ppm single element solutions as the NbOF_5^{-2} chemically stable for years in 2–5% HNO₃/trace HF in an LDPE container.

Nb CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HF/HNO₃); Oxide – very resistant to all acids including HF (fusion with K₂S₂O₇, KOH, or Na₂CO₃); Organic Matrices (dry ash at 450°C in Pt⁰ and dissolve by fusing with Na₂CO₃ or K₂S₂O₇).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 309.418 nm	0.04/0.002 µg/mL	1	ion	
ICP-OES 269.706 nm	0.07/0.002 µg/mL	1	ion	Th, Co
ICP-OES 295.088 nm	0.08/0.001 µg/mL	1	ion	Hf, U
ICP-MS 93 amu	1 ppt	n/a	M+	⁷⁷ Se ¹⁶ O, ⁷⁶ Se ¹⁷ O, [186X ⁺² (where X = W, Os)]

*ICP-OES D.L.'s are given as radial/axial view

Molybdenum

LOCATION: Group 6, Period 5

ATOMIC WEIGHT: 95.96

COORDINATION NUMBER: 6, 7, 8, 9 (coordination numbers 7, 8, and 9 are observed less frequently)

CHEMICAL FORM IN SOLUTION: MoO_4^{-2} (chem. form as received)

42	95.96
4639	1.3
2617	
Mo	
[Kr]4d ⁵ 5s	2,3,4,5,6
10.2	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Mo is received in a NH_4OH matrix, giving the operator the option of using HCl or HF to stabilize acidic solutions. The MoO_4^{-2} is soluble in concentrated HCl, MoOCl_5^{-2} , dilute HF/HNO_3 , MoOF_5^{-2} , and basic media MoO_4^{-2} . Stable at ppm levels with some metals, provided it is fluorinated. Do not mix with Alkaline or Rare Earths when HF is present. Stable with most inorganic anions, provided it is in the MoO_4^{-2} chemical form.

STABILITY: 2–100 ppb levels stable (alone or mixed with all other metals that are at comparable levels) as the MoOCl_5^{-2} for months in 1% HNO_3 /LDPE container. 1–10,000 ppm single element solutions as the MoO_4^{-2} chemically stable for years in 1% NH_4OH in a LDPE container.

Mo CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HF/HNO_3 or hot dilute HCl); Oxide (soluble in HF or NH_4OH); Organic Matrices (dry ash at 450°C in Pt^0 and dissolve oxide with HF or HCl).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 202.030 nm	0.008/0.0002 $\mu\text{g}/\text{mL}$	1	ion	Os, Hf
ICP-OES 203.844 nm	0.012/0.002 $\mu\text{g}/\text{mL}$	1	ion	
ICP-OES 204.598 nm	0.012/0.001 $\mu\text{g}/\text{mL}$	1	ion	Ir, Ta
ICP-MS 95 amu	3 ppt	n/a	M+	$^{40}\text{Ar}^{39}\text{K}^{16}\text{O}$, $^{79}\text{Br}^{16}\text{O}$, $^{190}\text{Os}^{2+}$, $^{190}\text{Pt}^{2+}$

*ICP-OES D.L.'s are given as radial/axial view

Ruthenium

LOCATION: Group 8, Period 5

ATOMIC WEIGHT: 101.07

COORDINATION NUMBER: 4, 5, 6, 8 (*coordination numbers 4, 5, and 8 are observed less frequently*)

CHEMICAL FORM IN SOLUTION: $[\text{RuCl}_6]^{2-}$

44	101.07
4150	1.4
2250	
Ru	
[Kr]4d ⁷ 5s	
12.4	2,3,4,6,8

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl. Stable with most metals and inorganic anions as the $[\text{RuCl}_6]^{2-}$ in dilute acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 10% HCl/LDPE container.

Ru CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (fuse with KOH/ KNO_3 in a Ag^0 crucible); Oxides (fuse with KOH/ KNO_3 in a Ag^0 crucible); Ores (see Oxides); Alloys (see Oxides). Organics (the RuO_4 is volatile and acidic oxidizing preparations should be used with caution. The preferred approach is the KOH/ KNO_3 fusion and dissolution of the fuseate in HCl).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 240.272 nm	0.03/0.002 $\mu\text{g/mL}$	1	ion	Fe
ICP-MS 101 amu	3 ppt	n/a	M+	⁴⁰ Ar ⁶¹ Ni, ⁶⁴ Ni ³⁷ Cl, ⁸⁵ Rb ¹⁶ O, ²⁰² Hg ²⁺
*ICP-OES D.L.'s are given as radial/axial view				

Rhodium

LOCATION: Group 6, Period 5

ATOMIC WEIGHT: 102.91

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: RhCl_6^{-3}

45
3697
1963
102.91
1.5

Rh

[Kr]4d⁸5s
12.4
2,3,4

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , H_2SO_4 and HF aqueous matrices. May cause AgCl precipitation when mixed with Ag^+ . Stable with all other metals.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 10% HCl/LDPE container.

Rh CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (elevated temp. with aqua regia or $\text{HCl}/\text{Cl}_2(\text{gas})$); Ores ($\text{HF}/\text{H}_2\text{SO}_4$ digestion followed by aqua regia digestion); Platinum scrap (aqua regia digestion).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 233.477 nm	0.04/0.004 $\mu\text{g}/\text{mL}$	1	ion	<i>Ni</i> , Sn, Mo, Nb, Ta
ICP-OES 249.077 nm	0.06/0.006 $\mu\text{g}/\text{mL}$	1	ion	Ta, Co, <i>Fe</i> , W, Cr, Os
ICP-OES 343.489 nm	0.06/0.006 $\mu\text{g}/\text{mL}$	1	atom	Mo, Th, Ce
ICP-MS 103 amu, monoisotopic	1 ppt	n/a	M+	⁴⁰ Ar ⁶³ Cu, ⁸⁷ Rb ¹⁶ O, ⁸⁷ Sr ¹⁶ O, ²⁰⁶ Pb ⁺²

*ICP-OES D.L.'s are given as radial/axial view

Palladium

LOCATION: Group 10, Period 5

ATOMIC WEIGHT: 106.42

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Pd}(\text{H}_2\text{O})_6^{2+}$

46	106.42
2964	1.4
1552	
Pd	
[Kr]4d ¹⁰	2,4
12.0	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl , HNO_3 , H_2SO_4 , HF , and H_3PO_4 . Avoid basic media. Stable with most metals and inorganic anions in acidic media. Avoid contact with water soluble organics such as aldehydes since Pd^{2+} is easily reduced.

STABILITY: 2–100 ppb levels. 2 ppb Pd is stable for 1 day in 1% HNO_3 /LDPE container. 10 ppb is stable for 3 days in 1% HNO_3 /LDPE container. 100 ppb is stable for ≥ 5 months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1-5% HNO_3 /LDPE container.

Pd CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3 or aqua regia); Oxides (soluble in HCl); Ores (dissolve in HCl/HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 340.458 nm	0.04/0.003 $\mu\text{g}/\text{mL}$	1	atom	Ce, Th, Zr
ICP-OES 363.470 nm	0.05/0.007 $\mu\text{g}/\text{mL}$	1	atom	
ICP-OES 229.651 nm	0.07/0.004 $\mu\text{g}/\text{mL}$	1	ion	Co
ICP-MS 105 amu	2 ppt	n/a	M+	$^{40}\text{Ar}^{65}\text{Cu}$, $^{89}\text{Y}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Silver

LOCATION: Group 11, Period 5

ATOMIC WEIGHT: 107.8682

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Ag}(\text{H}_2\text{O})_6^+$

47	107.868
2163	1.4
961	
Ag	
[Kr]4d ¹⁰ 5s	1
10.5	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HNO_3 and HF. Avoid basic media. Ag forms more insoluble salts than any other metal. It also is subject to photochemical reduction to the metal in HCl media although 10 $\mu\text{g}/\text{mL}$ solutions in 10% HCl $[\text{AgCl}_x^{1-x}]$ are commonly used in the analytical laboratory. The most common solubility problems exist with arsenate, arsenite, bromide, chloride, iodide, carbonate, chromate, cyanide, iodate, oxalate, oxide, sulfate, sulfide, tartrate, and thiocyanate in aqueous media. The addition of nitric acid renders many of these salts soluble.

STABILITY: 2–100 ppb levels stable for 75+ days when mixed with equivalent levels of all other elements including the precious metals (where chloride is present) when in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Ag CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); Oxides (soluble in HNO_3); Ores (digestion with conc. HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 328.068 nm	0.007/0.0007 $\mu\text{g}/\text{mL}$	1	atom	Ce, Rh, V
ICP-OES 338.289 nm	0.013/0.001 $\mu\text{g}/\text{mL}$	1	atom	Ce, Cr, Th
ICP-OES 243.779 nm	0.12/0.01 $\mu\text{g}/\text{mL}$	1	ion	Mn, Th, Ni, Rh
ICP-MS 107 amu	1 ppt	n/a	M+	$^{91}\text{Zr}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Cadmium

LOCATION: Group 12, Period 5

ATOMIC WEIGHT: 112.41

COORDINATION NUMBER: 4

CHEMICAL FORM IN SOLUTION: $\text{Cd}_2(\text{OH})(\text{aq})^{3+}$ and $\text{Cd}(\text{OH})(\text{aq})^{1+}$

48
112.41

767
1.5

321.18

Cd

[Kr]4d¹⁰5s²
2

8.65

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl, HNO_3 , H_2SO_4 , and HF. Avoid basic media forming insoluble carbonate and hydroxide. Stable with most metals and inorganic anions in acidic media. The sulfide, carbonate, oxalate, phosphate, and cyanide are insoluble in water and soluble in HCl, HNO_3 , and NH_4OH . The chloride, bromide, and iodide are soluble in water. CdI_2 is one of the few iodides soluble in ethanol. All compounds of Cd are soluble in excess NaI, due to the formation of the complex ion, CdI_4^{2-} .

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5 % HNO_3 /LDPE container.

Cd CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); Oxides (soluble in HCl or HNO_3); Ores (dissolve in HCl/ HNO_3 then take to fumes with H_2SO_4 . The silica and lead sulfate are filtered off after the addition of water); Organic based (dry ash at 450°C and dissolve ash in HCl), (sulfuric/peroxide acid digestion).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 214.438 nm	0.003/0.0003 $\mu\text{g/mL}$	1	ion	Pt, Ir
ICP-OES 228.802 nm	0.003/0.0003 $\mu\text{g/mL}$	1	atom	Co, Ir, As, Pt
ICP-OES 226.502 nm	0.003/0.0003 $\mu\text{g/mL}$	1	ion	Ir
ICP-MS 111 amu	11 ppt	n/a	M+	$^{95}\text{Mo}^{16}\text{O}$
*ICP-OES D.L.'s are given as radial/axial view				

Indium

LOCATION: Group 13, Period 5

ATOMIC WEIGHT: 114.82

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{In}(\text{H}_2\text{O})_6^{+3}$

49	114.82
2073	1.5
156.76	
In	
[Kr]4d ¹⁰ 5s ² 5p	
731	3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , and H_2SO_4 . Avoid neutral and basic media. Stable with most metals and inorganic anions. The oxalate, sulfide, carbonate, hydroxide, and phosphate are insoluble in water.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

In CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (is best dissolved in HCl/ HNO_3); Oxide (soluble in mineral acids); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (sulfuric/peroxide digestion or dry ash and dissolution in dilute HCl).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 158.583 nm	0.05/0.002 $\mu\text{g/mL}$	1	ion	
ICP-OES 230.606 nm	0.1/0.03 $\mu\text{g/mL}$	1	ion	<i>Ni, Os</i>
ICP-OES 325.609 nm	nm 0.2/0.05 $\mu\text{g/mL}$	1	atom	<i>Ir, Re</i>
ICP-MS 115 amu	1 ppt	n/a	M+	¹¹⁵ Sn, ⁹⁹ Ru ¹⁶ O

*ICP-OES D.L.'s are given as radial/axial view

Tin

LOCATION: Group 14, Period 5

ATOMIC WEIGHT: 118.710

COORDINATION NUMBER: 4, 5, 6, 7, 8 (coordination numbers 4, 5, 7 and 8 are observed less frequently)

CHEMICAL FORM IN SOLUTION: $\text{Sn}(\text{OH})_x\text{F}_y^{2-}$

50	118.71
2603	1.7
232.06	
Sn	
[Kr]4d ¹⁰ 5s ² 5p ²	
7.31	2,4

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and dilute HF/HNO₃. Avoid neutral to basic media. Unstable at ppm levels with metals that would pull F⁻ away (i.e. do not mix with Alkaline or Rare Earths or high levels of transition elements unless they are fluorinated). Stable with most inorganic anions, provided it is in the chemical form shown above.

STABILITY: 2–100 ppb levels stable (alone or mixed with all other metals that are at comparable levels) as the $\text{Sn}(\text{OH})_x\text{F}_y^{2-}$ for 1 year in 1% HNO₃/LDPE container. 1–10,000 ppm single element solutions as the $\text{Sn}(\text{OH})_x\text{F}_y^{2-}$ are chemically stable for years in 2–5% HNO₃/trace HF in a LDPE container.

Sn CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HF/HNO₃ or HCl); Oxides – SnO (soluble in HCl), SnO₂ – very resistant to all acids including HF (fusion with equal parts of Na₂CO₃ and S is soluble in water or dilute acids as the thiostannate); Alloys (treat first 0.1 g with 10 mL conc. H₂SO₄ to boiling until the alloy disintegrates and nearly all of the sulfuric acid is expelled. Then add 100 mL O₂ free water and 50 mL of conc. HCl or transfer to a plastic container and add 1 mL HF, in either case, warming gently to bring about solution); Organic Matrices (volatility and precipitation of the insoluble stannic oxide are problems – because these preparations are prone to error, we recommend you contact our technical staff at info@inorganicventures.com or (800) 669-6799 and we'll provide you with the necessary data for your specific sample type).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 189.989 nm	0.03/0.003 µg/mL	1	ion	
ICP-OES 242.949 nm	0.1/0.01 µg/mL	1	atom	W, Mo, Rh, Ta, Co
ICP-MS 120 amu	5 ppt	n/a	M+	¹²⁰ Te, ¹⁰⁴ Ru ¹⁶ O, ¹⁰⁴ Pd ¹⁶ O

*ICP-OES D.L.'s are given as radial/axial view

Antimony

LOCATION: Group 15, Period 5

ATOMIC WEIGHT: 121.760

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Sb}(\text{O})\text{C}_4\text{H}_4\text{O}_6^{-1}$

51
1587
6317
121.76
1.8

Sb

[Kr]4d¹⁰5s²5p³
6.6
±3,5

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in conc. HCl, dilute or conc. HF. Stable in dilute HNO_3 as the fluoride or tartrate complex. Avoid basic media. Stable with most metals and inorganic anions in acidic media as the tartrate provided the acidity is not too high or the acid is oxidizing causing loss of the stabilizing tartrate ion. The fluoride complex of antimony is stable in strong acid but you should only mix with other metals that are fluorinated.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–2% HNO_3 /LDPE container.

Sb CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal and alloys (soluble in H_2O /HF/ HNO_3 mixture); Oxides (soluble in HCl and tartaric acid or H_2O /HF/ HNO_3 mixtures); Ores (fusion with Na_2CO_3 in Pt^0 followed by dissolving the fuseate in a H_2O /HF/ HNO_3 mixture); Organic based (sulfuric acid/hydrogen peroxide digestion).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 206.833 nm	0.03/0.003 $\mu\text{g}/\text{mL}$	1	atom	<i>Ta, Cr, Ge, Hf</i>
ICP-OES 217.581 nm	0.05/0.005 $\mu\text{g}/\text{mL}$	1	atom	<i>Nb, W, Re, Fe</i>
ICP-OES 231.147 nm	0.06/0.006 $\mu\text{g}/\text{mL}$	1	atom	<i>Ni, Co, Pt</i>
ICP-MS 121 amu	5 ppt	n/a	M+	¹⁰⁵ Pd ¹⁶ O, ⁸⁹ Y ¹⁶ O ₂

*ICP-OES D.L.'s are given as radial/axial view

Tellurium

LOCATION: Group 16, Period 5

ATOMIC WEIGHT: 127.60

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: H_2TeO_3 (HNO_3 matrix), TeCl_6^{2-} (HCl matrix)

52	127.60
988	2.0
449.65	
Te	
[Kr]4d ¹⁰ 5s ² 5p ⁴	
6.24	-2,4,6

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl , HNO_3 , H_3PO_4 , H_2SO_4 and HF aqueous matrices and water. It is stable with most inorganic anions and cations. Avoid mixing HCl matrices with elements forming insoluble chlorides such as Ag^+ . When fluorinated and/or under acidic conditions precipitation is typically not a problem at moderate to low concentrations.

STABILITY: 2–100 ppb levels stable for months alone or mixed with other elements at equivalent levels in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Te CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in solutions of alkali hydroxides or a 1:1:1 mixture of H_2O , H_2SO_4 , HNO_3); Oxides (TeO_2 is soluble in HCl and the alkali hydroxides. TeO_3 is soluble in hot concentrated solutions of the alkali hydroxides.); Minerals and alloys (acid digestion with HNO_3 or HNO_3 / HF); Organic Matrices (Vegetable Matter — dry ash 100 g of the well-ground and mixed vegetation into a concentrated solution of 25 g of magnesium nitrate and magnesium oxide. Dry, ignite and muffle until the ash is a uniform gray color).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 170.000 nm	0.04/0.004 $\mu\text{g/mL}$	1	atom	Sn
ICP-OES 214.281 nm	0.04/0.004 $\mu\text{g/mL}$	1	atom	Ta, Re, V
ICP-OE 225.902 nm	0.20/0.02 $\mu\text{g/mL}$	1	atom	Ir, Os W, Ga, Ru, Ta
ICP-MS 130 amu	20 ppt	n/a	M+	$^{114}\text{Cd}^{16}\text{O}$, $^{114}\text{Sn}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Cesium

LOCATION: Group 1, Period 6

ATOMIC WEIGHT: 132.9054

COORDINATION NUMBER: (6) (*coordination number in parentheses is assumed, not certain*)

CHEMICAL FORM IN SOLUTION: Cs+(aq)

55	132.905
671	0.9
28.55	
Cs	
[Xe]6s	1
1.873	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in dilute HCl, HNO_3 , H_2SO_4 and HF aqueous matrices. Stable with most metals and inorganic anions. Forms insoluble $\text{Cs}_2[\text{PtCl}_6]$.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 0.1% HNO_3 /LDPE container. Will crystallize out of higher (~ 5%) levels of HNO_3 at > 1000 $\mu\text{g/mL}$.

Cs CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (dissolves very rapidly in water); Ores (sodium carbonate fusion in Pt^0 followed by HCl dissolution – blank levels of Cs in sodium carbonate critical); Organic Matrices (sulfuric/peroxide digestion or nitric/sulfuric/perchloric acid decomposition).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 455.531 nm	100/2 $\mu\text{g/mL}$	1	atom	Cr, U, Ce, Ti
ICP-MS 133 amu	1.7 ppt	n/a	M+	$^{117}\text{Sn}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Barium

LOCATION: Group 2, Period 6

ATOMIC WEIGHT: 137.33

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Ba}(\text{H}_2\text{O})_6^{+2}$

56	137.33
1898	1.0
729	
Ba	
[Xe]6s ²	2
3.5	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid H_2SO_4 , HF, and neutral to basic media. Stable with most metals and inorganic anions forming insoluble silicate, carbonate, hydroxide, oxide, fluoride, sulfate, oxalate, chromate, arsenate, iodate, molybdate, sulfite and tungstate in neutral aqueous media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Ba CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (is best dissolved in diluted HNO_3); Ores (carbonate fusion in Pt^0 followed by HCl dissolution. If sulfate is present dissolve the fuseate using HCl/tartaric acid to prevent BaSO_4 precipitate); Organic Matrices (dry ash and dissolve in dilute HCl. Do not heat when dissolving to avoid precipitation of SiO_2).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 455.403 nm	0.002/0.0001 $\mu\text{g/mL}$	1	ion	Zr, U
ICP-OES 233.527 nm	0.004/0.0003 $\mu\text{g/mL}$	1	ion	
ICP-OES 230.424 nm	0.004/0.0005 $\mu\text{g/mL}$	1	ion	Mo, Ir, Co
ICP-MS 138 amu	1 ppt	n/a	M+	$^{122}\text{Sn}^{16}\text{O}$, $^{122}\text{Te}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Lanthanum

LOCATION: Group 3, Period 6 (lanthanoid)

ATOMIC WEIGHT: 138.9055

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{La}(\text{OH})_y(\text{H}_2\text{O})_x^{+3-y}$

57	138.906
3457	1.1
920	
La	
[Xe]5d6s ²	
6.145	3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid HF, H_3PO_4 , H_2SO_4 and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride and sparingly soluble sulfates (La–Eu exhibit low sulfate solubility). Avoid mixing with elements / solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

La CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 333.749 nm	0.01/0.001 $\mu\text{g}/\text{mL}$	1	ion	
ICP-OES 408.672 nm	0.01/0.001 $\mu\text{g}/\text{mL}$	1	ion	Th
ICP-OES 412.323 nm	0.01/0.001 $\mu\text{g}/\text{mL}$	1	ion	Ce, Th
ICP-MS 139 amu	1 ppt	n/a	M+	$^{123}\text{Sb}^{16}\text{O}$, $^{123}\text{Te}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Hafnium

LOCATION: Group 4, Period 6

ATOMIC WEIGHT: 178.49

COORDINATION NUMBER: 6, 7, 8 (*coordination numbers 7 & 8 are observed less frequently*)

CHEMICAL FORM IN SOLUTION: $\text{Hf}(\text{F})_6^{-2}$

72	178.49
4603	1.2
2227	
Hf	
[Xe]4f ¹⁴ 5d ² 6s ²	
13.3	4

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in concentrated HCl, HF, H_2SO_4 (very hot), and HNO_3 . Avoid H_3PO_4 and neutral to basic media. Unstable at ppm levels with metals that would pull F⁻ away (i.e. do not mix with Alkaline or Rare Earths or high levels of transition elements unless they are fluorinated). Stable with most inorganic anions but precipitation with phosphate, oxalate, and tartrate with a tendency to hydrolyze forming the hydrated oxide in all dilute acids except HF.

STABILITY: 2–100 ppb levels stable alone or mixed with all other metals that are at comparable levels as the $\text{Hf}(\text{F})_6^{-2} + \text{Hf}(\text{OH})_4\text{F}_2^{-2}$ for months in 1% HNO_3 /LDPE container. 1–10,000 ppm single element solutions as the $\text{Hf}(\text{F})_6^{-2}$ chemically stable for years in 2-5% HNO_3 /trace HF in an LDPE container.

HF CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in $\text{H}_2\text{O}/\text{HF}/\text{HNO}_3$); Oxide — unlike TiO_2 the HfO_2 is best fused in one of the following ways (Na_2O_2 in Ni^0 , Na_2CO_3 in Pt^0 or Borax in Pt^0); Organic Matrices (dry ash at 450°C in Pt^0 and dissolve by fusing with Na_2CO_3 and dissolving in $\text{HF}/\text{HNO}_3/\text{H}_2\text{O}$).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 277.336 nm	0.02/0.002 µg/mL	1	ion	Nb, Cr, U
ICP-OES 273.876 nm	0.02/0.002 µg/mL	1	ion	U, Mo
ICP-OES 264.141 nm	0.02/0.002 µg/mL	1	ion	Ba, Th, U
ICP-MS 177 amu	4 ppt	n/a	M+	¹⁶¹ Dy ¹⁶ O**

*ICP-OES D.L.'s are given as radial/axial view

**Fewer potential interferences on the 177 vs 180 mass

Tantalum

LOCATION: Group 5, Period 6

ATOMIC WEIGHT: 180.9479

COORDINATION NUMBER: 6, 7, 8 (*coordination numbers 7 & 8 are observed less frequently*)

CHEMICAL FORM IN SOLUTION: TaOF_6^{-3}

73	180.948
5458	1.3
3014	
Ta	
[Xe]4f ¹⁴ 5d ³ 6s ²	
16.6	5

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in concentrated HCl and dilute HF/HNO₃. Avoid neutral to basic media. Unstable at ppm levels with metals that would pull F⁻ away (i.e. do not mix with Alkaline or Rare Earths or high levels of transition elements unless they are fluorinated). Stable with most inorganic anions, provided it is in the chemical form shown above.

STABILITY: 2–100 ppb levels stable (alone or mixed with all other metals that are at comparable levels) as the TaOF_6^{-3} , for 2 months at the 2–10 ppb level in 1% HNO₃/LDPE container and for 5 months at the 100 ppb level under same conditions. 1–10,000 ppm single element solutions as the TaOF_6^{-3} are chemically stable for years in 2–5% HNO₃/trace HF in an LDPE container.

Ta CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HF/HNO₃); Oxide — very resistant to all acids including HF (fusion with K₂S₂O₇, KOH, or Na₂CO₃); Organic Matrices (dry ash at 450°C in Pt⁰ and dissolve by fusing with Na₂CO₃ or K₂S₂O₇).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 226.230 nm	0.03/0.01 µg/mL	1	ion	Sb, Nb
ICP-OES 240.063 nm	0.03/0.004 µg/mL	1	ion	Hf, Fe, Bi
ICP-OES 268.517 nm	0.03/0.005 µg/mL	1	ion	Cr, Ru, Hf, W
ICP-MS 181 amu	2 ppt	n/a	M+	¹⁶⁵ Ho ¹⁶ O

*ICP-OES D.L.'s are given as radial/axial view

Tungsten

LOCATION: Group 6, Period 6

ATOMIC WEIGHT: 183.84

COORDINATION NUMBER: 6, 7, 8, 9 (*coordination numbers are 7, 8 and 9 are observed less frequently*)

CHEMICAL FORM IN SOLUTION: WOF_5^{-2} (chem. form as received)

74
5555
3407
183.84
1.4

W

[Xe]4f¹⁴5d⁴6s²
19.3
2,3,4,5,6

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: W is very readily hydrolyzed requiring 0.1 to 1% HF for stable acidic solutions. The $[\text{WOF}_5]^{-2}$ is soluble in % levels of HCl and HNO_3 , provided it is in the $[\text{WOF}_5]^{-2}$ form. Stable at ppm levels with some metals provided it is fluorinated. Do not mix with Alkaline or Rare Earths. W is best to be mixed only with other fluorinated metals (Ti, Zr, Hf, Nb, Ta, Mo, Si, Sn, Ge). Look for yellow WO_3 precipitate if mixed with other transitions at higher levels indicating instability. The yellow WO_3 will form over a period of weeks even in trace HF, *therefore HF levels of W multi-element blends should be ~ 1%.*

STABILITY: 2–100 ppb levels stable (alone or mixed with all other metals that are at comparable levels) as the $[\text{WOF}_5]^{-2}$ for months in 1% HNO_3 /LDPE container. 1–10,000 ppm single element solutions as the $[\text{WOF}_5]^{-2}$ chemically stable for years in 1% HF in a LDPE container.

W CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HF/ HNO_3); Oxide (soluble in HF or NH_4OH); Organic Matrices (dry ash at 450°C in Pt^0 and dissolve oxide with HF).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 207.911 nm	0.03/0.001 $\mu\text{g/mL}$	1	ion	Ru, In
ICP-OES 224.875 nm	0.05/0.005 $\mu\text{g/mL}$	1	ion	Co, Rh, Ag
ICP-OES 209.475 nm	0.05/0.008 $\mu\text{g/mL}$	1	ion	Mo
ICP-MS 182 amu	5 ppt	n/a	M+	$^{166}\text{Er}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Rhenium

LOCATION: Group 8, Period 6

ATOMIC WEIGHT: 186.207

COORDINATION NUMBER: 4, 6, 7, 8, 9 (coordination numbers 4, 7, 8 and 9 are observed less frequently)

CHEMICAL FORM IN SOLUTION: ReO_4^{1-}

75
5596
3180

186.207
1.5

Re

[Xe]4f¹⁴5d⁵6s²
21.0
-1,2,4,6,7

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl , HNO_3 , H_2SO_4 , HF , and H_3PO_4 . Stable with most metals and inorganic anions in acidic media. Mixing higher levels of ReO_4^{1-} with Ag^+ , Hg_2^{2+} , K^+ , NH_4^+ , Cs^+ , Rb^+ , or Tl^+ will give the corresponding salt (solubilities are 1–12 g/L).

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 1–5% HNO_3 /LDPE container.

Re CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); Oxides/Ores (fuse in Pt^0 with Na_2CO_3). Organic Matrix (all modes of acid attack invite the danger of loss of some volatile perrhenic acid. The use of a reflux condenser should be considered when a wet acid digestion is used such as nitric/perchloric or sulfuric/peroxide digestions. The preferred approach is to ash the sample in Pt^0 mixed with Na_2CO_3 starting the ash at 450°C and then increasing the temperature, if necessary, to 900°C to effect a fusion of accompanying aluminosilicates, etc.).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 221.426 nm	0.006/0.0006 $\mu\text{g/mL}$	1	ion	Fe, Os, Mo, Ta
ICP-OES 227.525 nm	0.006/0.0006 $\mu\text{g/mL}$	1	ion	Ru, Co, Ca
ICP-MS 187 amu	2 ppt	n/a	M+	$^{171}\text{Yb}^{16}\text{O}$, ^{187}Os
*ICP-OES D.L.'s are given as radial/axial view				

Osmium

LOCATION: Group 8, Period 6

ATOMIC WEIGHT: 190.23

COORDINATION NUMBER: 4, 5, 6, 8 (*coordination numbers are 4, 5 and 8 are observed less frequently*)

CHEMICAL FORM IN SOLUTION: OsCl_6^{2-}

76	190.23
5012	1.5
3027	
Os	
[Xe]4f ¹⁴ 5d ⁶ 6s ²	
22.6	2,3,4,6,8

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl. Stable with most metals and inorganic anions as the $[\text{OsCl}_6]^{2-}$ in dilute HCl media. **DO NOT EXPOSE TO NITRIC ACID – FORMATION OF THE VERY VOLATILE AND TOXIC OsO_4 WILL RESULT.** Any oxidizing condition must be avoided.

STABILITY: 2–100 ppb levels are *NOT* stable in 1% HNO_3 /LDPE container. The stability of HCl solutions at ppb levels has not been determined by our laboratory. 1–10,000 ppm solutions are presumed chemically stable for years in 10% HCl/LDPE container, stability studies have not been performed.

Os CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal powder – Fuse with NaOH/ NaNO_3 in quartz crucibles at approximately 400C. See, *Tracing an Osmium Solution Standard to the International System of Units (SI)*, Madeline Gozzi, et al., Anal. Chem. 2021, 93, 15642–15650; Ores & Oxides – Fuse with KOH/ KNO_3 in a Ag₂O crucible and dissolve in water being sure to avoid addition of any acid; Organics – the OsO_4 is volatile and acidic oxidizing preparations should be used with caution.

NOTE: The presence of the OsO_4 will give false high results due to its enhanced nebulization efficiency (volatility). *Only dilutions in HCl should be made. Adding 0.5% Hydroxylamine Hydrochloride to sample preparations will help with ICP washout. The use of HNO_3 should be strictly avoided.* Preparations from caustic nitrate fusions should be diluted in water.

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 225.585 nm	0.0004 $\mu\text{g/mL}$	1	ion	Fe, Ta, Ge, Ir, Cr
ICP-MS 192 amu	1 ppt	n/a	M+	$^{176}\text{Yb}^{16}\text{O}$, $^{176}\text{Lu}^{16}\text{O}$, $^{176}\text{Hf}^{16}\text{O}$, ^{192}Pt

*ICP-OES D.L.'s are given as radial/axial view

Iridium

LOCATION: Group 9, Period 6

ATOMIC WEIGHT: 192.22

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: IrCl_6^{-2}

77	192.22
4428	1.6
2443	
Ir	
[Xe]4f ¹⁴ 5d ⁷ 6s ²	
22.4	2,3,4,6

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , H_2SO_4 , and HF aqueous matrices. May cause AgCl precipitation when mixed with Ag^+ . Stable with all other metals.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 10% HCl/LDPE container.

Ir CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (Fusion) — Heat a mixture of Ir metal powder and CaCO_3 at 1200°C for 3 hours. The Ir metal must be a very fine mesh powder to react completely with the CaCO_3 . Convert the Ca_4IrO_6 to IrCl_4 and bring into solution by heating almost to boiling for a least 1 hour in 80% HCl; Metal — Heat with aqua regia or HCl/ Cl_2 (gas); Ores — HF/ H_2SO_4 digestion followed by aqua regia digestion; Platinum scrap — Aqua regia digestion.

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 224.268 nm	0.03 $\mu\text{g}/\text{mL}$	1	ion	<i>Cu, Nb, Hf</i>
ICP-OES 212.681 nm	0.03 $\mu\text{g}/\text{mL}$	1	ion	<i>Ta, Yb, Au, V</i>
ICP-OES 205.222 nm	0.06 $\mu\text{g}/\text{mL}$	1	atom	<i>Fe</i>
ICP-MS 191 amu	2 ppt	n/a	M+	¹⁷⁵ <i>Lu</i> ¹⁶ O

*ICP-OES D.L.'s are given as radial/axial view

Platinum

LOCATION: Group 10, Period 6

ATOMIC WEIGHT: 195.078

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Pt}(\text{Cl})_6^{2-}$

78
3827
1772
195.08
1.4

Pt

[Xe]4f¹⁴5d⁶6s²
22.6
2,3,4,6,8

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl and HNO_3 , as the chloride complex. Avoid basic media. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–10 ppb Pt is stable for 2 months in 1% HNO_3 /LDPE container. 100 ppb is stable for 5 months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 10 % HCl/LDPE container.

Pt CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (aqua regia); Oxides (soluble in HCl); Ores (dissolve in HCl/ HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 214.423 nm	0.03/0.003 $\mu\text{g/mL}$	1	ion	W, As, Ir, Cd
ICP-OES 203.646 nm	0.06/0.006 $\mu\text{g/mL}$	1	ion	Co, Hf
ICP-MS 195 amu	5 ppt	n/a	M+	¹⁷⁹ Hf ¹⁶ O

*ICP-OES D.L.'s are given as radial/axial view

Gold

LOCATION: Group 11, Period 6

ATOMIC WEIGHT: 196.9665

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Au}(\text{Cl})_6^{3-}$

79	196.967
2857	1.4
1064.58	
Au	
[Xe]4f ¹⁴ 5d ¹⁰ 6s	
18.9	1,3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl, and HNO_3 , as the chloride complex. Avoid basic media. Stable with most metals and inorganic anions in acidic media.

STABILITY: 2–100 ppb levels. 2-10 ppb Au is stable for ≤ 1 day maximum in 1% HNO_3 /LDPE container. 100 ppb is stable for ≤ 2 days maximum in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 10% HCl/LDPE container.

Au CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (aqua regia); Oxides (soluble in HCl); Ores (dissolve in HCl/ HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 242.795 nm	0.02/0.003 $\mu\text{g/mL}$	1	atom	<i>Mn, Os, Th, Ta, Pt</i> Co, F
ICP-OES 267.595 nm	0.03/0.003 $\mu\text{g/mL}$	1	atom	<i>Nb, Ta, U, Cr, Th,</i> Rh, Ru
ICP-OES 208.209 nm	0.04/0.01 $\mu\text{g/mL}$	1	ion	Ir, Re
ICP-MS 197 amu	5 ppt	n/a	M+	¹⁸¹ Ta ¹⁶ O

*ICP-OES D.L.'s are given as radial/axial view

Mercury

LOCATION: Group 12, Period 6

ATOMIC WEIGHT: 200.59

COORDINATION NUMBER: 4

CHEMICAL FORM IN SOLUTION: $\text{Hg}(\text{OH})(\text{aq})^{1+}$

80	200.59
357	1.5
-38.72	
Hg	
[Xe]4f ¹⁴ 5d ⁶ 6s ²	
13.5	1,2

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HNO_3 . Avoid basic media forming insoluble carbonate. The sulfide, basic carbonate, oxalate, phosphate, arsenite, arsenate, and iodide are insoluble in water.

STABILITY: 2–100 ppb levels — stable in 10% HNO_3 packaged in borosilicate glass; *NOT* stable in 1% HNO_3 /LDPE container. 1–100 ppm levels stable in 7% HNO_3 packaged in borosilicate glass. 1000–10,000 ppm solutions are chemically stable for years in 5–10% HNO_3 /LDPE container.

Hg CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); HgO (soluble in HNO_3); Ores and Organic based (our documentation has more references to the preparation of Hg containing samples than any other element — because these preparations are prone to error, we recommend you contact our technical staff at info@inorganicventures.com or (800) 669-6799 and we'll provide you with the necessary data for your specific sample type).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 184.950 nm	0.03/0.005 $\mu\text{g/mL}$	1	atom	
ICP-OES 194.227 nm	0.03/0.005 $\mu\text{g/mL}$	1	ion	V
ICP-OES 253.652 nm	0.1 /0.03 $\mu\text{g/mL}$	1	atom	Ta, Co, Th ,Rh , Fe, U
ICP-MS 202 amu	9 ppt	n/a	M+	¹⁸⁶ W ¹⁶ O

*ICP-OES D.L.'s are given as radial/axial view

Thallium

LOCATION: Group 13, Period 6

ATOMIC WEIGHT: 204.383

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Tl}(\text{H}_2\text{O})_6^{1+}$

81	204.383
1473	1.4
304	
TI	
[Xe]4f ¹⁴ 5d ¹⁰ 6s ² 6p ^{1,3}	
11.85	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HNO_3 , and H_2SO_4 . Stable with most metals and inorganic anions. The sulfite, thiocyanate, and oxalate are moderately soluble; the phosphate and arsenite are slightly soluble and the sulfide is insoluble.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

TI CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (is best dissolved in HNO_3 which forms chiefly the Tl^{1+} ion); Oxide (the thallic oxide is readily soluble in water. The thallic oxide requires high levels of acid); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (sulfuric/peroxide digestion or dry ash and dissolution in HCl).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 190.864 nm	0.04/0.004 $\mu\text{g/mL}$	1	ion	V, Ti
ICP-OES 276.787 nm	0.1/0.01 $\mu\text{g/mL}$	1	atom	Ta, V, Fe, Cr
ICP-OES 351.924 nm	0.2/0.02 $\mu\text{g/mL}$	1	atom	Th, Ce, Zr
ICP-MS 205 amu	2 ppt	n/a	M+	$^{189}\text{Os}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Lead

LOCATION: Group 14, Period 6

ATOMIC WEIGHT: 207.2

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Pb}(\text{H}_2\text{O})_6^{+2}$

82
207.2

1750
1.6

Pb

[Xe]4f¹⁴5d¹⁰6s²6p²
13.5

2,4

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, HF, and HNO_3 . Avoid H_2SO_4 . Stable with most metals and inorganic anions forming insoluble carbonate, borate, sulfate, sulfite, sulfide, phosphate, oxalate, chromate, tannate, iodate, and cyanide in neutral aqueous media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Pb CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (best dissolved in 1:1 $\text{H}_2\text{O}/\text{HNO}_3$); Oxides (the many different Pb oxides are soluble in HNO_3 , with the exception of PbO_2 which is soluble in HCl or HF); Ores and Alloys (best attacked using 1:1 $\text{H}_2\text{O}/\text{HNO}_3$); Organic Matrices (dry ash and dissolve in dilute HCl. Do not heat when dissolving to avoid precipitation of SiO_2).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 168.215 nm	0.03/0.003 $\mu\text{g}/\text{mL}$	1	ion	Co
ICP-OES 220.353 nm	0.04/0.006 $\mu\text{g}/\text{mL}$	1	ion	Bi, Nb
ICP-OES 217.000 nm	0.09/0.03 $\mu\text{g}/\text{mL}$	1	atom	W, Ir, Hf, Sb, Th
ICP-MS 208 amu	5 ppt	n/a	M+	$^{192}\text{Pt}^{16}\text{O}$, $^{192}\text{Os}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Bismuth

LOCATION: Group 15, Period 6

ATOMIC WEIGHT: 208.9804

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Bi}(\text{O})(\text{H}_2\text{O})_x^{1+}$

83
208.980

1564
1.4

271.52

Bi

[Xe]4f¹⁴5d¹⁰6s²6p³
9.75

3,5

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Stable in HCl , HNO_3 , H_2SO_4 , and HF . Avoid basic media forming insoluble hydroxide. Stable with most metals and inorganic anions in acidic media. Many salts that are insoluble in water are soluble in HCl , HNO_3 and HF . The major problem with Bi^{3+} is its tendency to hydrolyze at higher concentrations or in dilute acid. Nitric acid solutions should be 5% to hold the Bi in solution in the 100 to 10000 $\mu\text{g}/\text{mL}$ concentration range.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 5–7% HNO_3 /LDPE container.

BI CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in HNO_3); Oxides (soluble in HNO_3); Alloys (dissolve in conc. 4:1 HCl/HNO_3 — heating may be required.); Organic based (dry ash at 450°C and dissolve ash in HNO_3 or acid digestion with conc. hot sulfuric acid adding hydrogen peroxide *carefully dropwise* until clear).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 223.061 nm	0.04/0.005 $\mu\text{g}/\text{mL}$	1	atom	Th, Ir, Ti Cu
ICP-OES 306.772 nm	0.08/0.01 $\mu\text{g}/\text{mL}$	1	atom	<i>Th</i> , U, Zr, Hf, Fe
ICP-OES 222.825 nm	0.1/0.02 $\mu\text{g}/\text{mL}$	1	atom	<i>Cr</i> , <i>Hf</i> , Ce, Os
ICP-MS 209 amu	2 ppt	n/a	M+	¹⁹³ Ir ¹⁶ O

*ICP-OES D.L.'s are given as radial/axial view

Cerium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 140.12

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Ce}(\text{OH})_y(\text{H}_2\text{O})_x^{+4-y}$

58	140.12
3426	1.1
798	
Ce	
[Xe]4f ² 6s ²	3,4
6.66	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid HF, H_3PO_4 , H_2SO_4 and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride and sparingly soluble sulfates (La–Eu exhibit low sulfate solubility). Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Ce CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 413.765 nm	0.05/0.0058 $\mu\text{g}/\text{mL}$	1	ion	Ce**
ICP-OES 418.660 nm	0.05/0.003 $\mu\text{g}/\text{mL}$	1	ion	Zr
ICP-OES 453.975 nm	0.06/0.0063 $\mu\text{g}/\text{mL}$	1	ion	
ICP-MS 140 amu	1 ppt	n/a	M+	$^{124}\text{Sn}^{16}\text{O}$, $^{124}\text{Te}^{16}\text{O}$
*ICP-OES D.L.'s are given as radial/axial view **413.747 line may effect Bkg. Corr.				

Praseodymium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 140.9077

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Pr}(\text{OH})_y(\text{H}_2\text{O})_x^{+3-y}$

59	140.908
3512	1.1
931	
Pr	
[Xe]4f ³ 6s ²	3,4
6.77	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid HF, H_3PO_4 , H_2SO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride and sparingly soluble sulfates (La–Eu exhibit low sulfate solubility). Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Pr CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 414.311 nm	0.04/0.004 $\mu\text{g}/\text{mL}$	1	ion	Ce
ICP-OES 417.939 nm	0.04/0.004 $\mu\text{g}/\text{mL}$	1	ion	Cr, Ce
ICP-OES 422.535 nm	0.04/0.004 $\mu\text{g}/\text{mL}$	1	ion	V, U
ICP-MS 141 amu	0.3 ppt	n/a	M+	$^{125}\text{Te}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Neodymium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 144.24

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Nd}(\text{OH})_y(\text{H}_2\text{O})_x^{+3-y}$

60	144.24
3068	1.1
1016	
Nd	
[Xe]4f ⁴ 6s ²	
7.00	3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid HF, H_3PO_4 , H_2SO_4 and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride and sparingly soluble sulfates (La–Eu exhibit low sulfate solubility). Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Nd CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 401.225 nm	0.05/0.002 $\mu\text{g}/\text{mL}$	1	ion	Ti, Cr
ICP-OES 430.358 nm	0.075/0.0014 $\mu\text{g}/\text{mL}$	1	ion	
ICP-OES 406.109 nm	0.1/0.002 $\mu\text{g}/\text{mL}$	1	ion	Ce
ICP-MS 146 amu	2 ppt	n/a	M+	$^{130}\text{Te}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Samarium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 150.36

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Sm}(\text{OH})_y(\text{H}_2\text{O})_x^{+3-y}$

62	150.36
1791	1.1
1072	
Sm	
[Xe]4f ⁶ 6s ²	2,3
7.52	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid HF, H_3PO_4 , H_2SO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride and sparingly soluble sulfates (La–Eu exhibit low sulfate solubility). Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Sm CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 359.260 nm	0.05/0.002 $\mu\text{g}/\text{mL}$	1	ion	W, Th
ICP-OES 442.434 nm	0.075/0.0014 $\mu\text{g}/\text{mL}$	1	ion	Ce, Ca
ICP-OES 428.079 nm	0.1/0.002 $\mu\text{g}/\text{mL}$	1	ion	Ce, Cr
ICP-MS 152 amu	2 ppt	n/a	M+	$^{136}\text{Ce}^{16}\text{O}$, $^{136}\text{Ba}^{16}\text{O}$, ^{152}Gd

*ICP-OES D.L.'s are given as radial/axial view

Europium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 151.96

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Eu}(\text{OH})_y(\text{H}_2\text{O})_x^{+3-y}$

63	151.96
1597	1.0
817	
Eu	
[Xe]4f ⁷ 6s ²	2,3
5.24	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid HF, H_3PO_4 , H_2SO_4 and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride and sparingly soluble sulfates (La–Eu exhibit low sulfate solubility). Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Eu CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 381.967 nm	0.003/0.0003 $\mu\text{g}/\text{mL}$	1	ion	Cr, V
ICP-OES 412.970 nm	0.004/0.0004 $\mu\text{g}/\text{mL}$	1	ion	Nb
ICP-OES 420.505 nm	0.004/0.0004 $\mu\text{g}/\text{mL}$	1	ion	Ce, V
ICP-MS 153 amu	1 ppt	n/a	M+	$^{137}\text{Ba}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Gadolinium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 157.25

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Gd}(\text{OH})_x(\text{H}_2\text{O})_y^{+3-x}$

64	157.25
3266	1.1
1312	
Gd	
[Xe]4f ⁷ 5d6s ²	
7.92	3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl , H_2SO_4 , and HNO_3 . Avoid HF , H_3PO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Gd CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 342.247 nm	0.014/0.001 $\mu\text{g}/\text{mL}$	1	ion	Th, U
ICP-OES 336.223 nm	0.02/0.0002 $\mu\text{g}/\text{mL}$	1	ion	Th, Ca
ICP-OES 335.047 nm	0.02/0.002 $\mu\text{g}/\text{mL}$	1	ion	Ce, Ca
ICP-MS 158 amu	2 ppt	n/a	M+	$^{142}\text{Ce}^{16}\text{O}$, $^{142}\text{Nd}^{16}\text{O}$, ^{158}Dy
*ICP-OES D.L.'s are given as radial/axial view				

Terbium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 158.925

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Tb}(\text{OH})_x(\text{H}_2\text{O})_y^{+3-x}$

65	158.925
3223	1.1
1357	
Tb	
[Xe]4f ⁹ 6s	3,4
8.23	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, H_2SO_4 , and HNO_3 . Avoid HF, H_3PO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Tb CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 350.917 nm	0.02/0.002 $\mu\text{g}/\text{mL}$	1	ion	V, Th, Ce, Zr
ICP-OES 367.635 nm	0.06/0.006 $\mu\text{g}/\text{mL}$	1	ion	Ta, Ce, Co, U
ICP-MS 159 amu	1 ppt	n/a	M+	$^{143}\text{Nd}^{16}\text{O}$, $^{127}\text{I}^{16}\text{O}_2$

*ICP-OES D.L.'s are given as radial/axial view

Dysprosium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 162.50

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Dy}(\text{OH})_x(\text{H}_2\text{O})_y^{+3-x}$

66	162.50
2562	1.1
1409	
Dy	
[Xe]4f ¹⁰ 6s ²	
8.55	3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, H₂SO₄, and HNO₃. Avoid HF, H₃PO₄, and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO₃/LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO₃/LDPE container.

Dy CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in H₂O/HNO₃); Ores (carbonate fusion in Pt⁰ followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 H₂O/HCl or HNO₃).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 340.780 nm	0.007/0.0007 µg/mL	1	ion	Hf, Th, U, Zr
ICP-OES 353.170 nm	0.013/0.001 µg/mL	1	ion	Ce, Th
ICP-MS 163 amu	3 ppt	n/a	M+	¹⁴⁷ Sm ¹⁶ O
*ICP-OES D.L.'s are given as radial/axial view				

Holmium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 164.930

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Ho}(\text{OH})_x(\text{H}_2\text{O})_y^{+3-x}$

67	164.93
2695	1.1
1470	
Ho	
[Xe]4f ¹¹ 6s ²	3
8.80	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, H_2SO_4 , and HNO_3 . Avoid HF, H_3PO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Ho CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 345.600 nm	0.006/0.0001 µg/mL	1	ion	U, Ti
ICP-OES 339.898 nm	0.02/0.002 µg/mL	1	ion	Ce, Re
ICP-MS 165 amu	1 ppt	n/a	M+	$^{149}\text{Sm}^{16}\text{O}$

*ICP-OES D.L.'s are given as radial/axial view

Erbium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 167.26

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Er}(\text{OH})_x(\text{H}_2\text{O})_y^{+3-x}$

68	167.26
2863	1.1
1522	
Er	
[Xe]4f ¹² 6s ²	
9.07	3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, H₂SO₄, and HNO₃. Avoid HF, H₃PO₄, and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO₃/LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO₃/LDPE container.

Er CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in H₂O/HNO₃); Ores (carbonate fusion in Pt⁰ followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 H₂O/HCl or HNO₃).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 337.271 nm	0.01/0.001 µg/mL	1	ion	Th, Ti
ICP-OES 349.910 nm	0.02/0.002 µg/mL	1	ion	Ru, Th, U
ICP-MS 166 amu	1 ppt	n/a	M+	¹⁵⁰ Sm ¹⁶ O, ¹⁵⁰ Nd ¹⁶ O
*ICP-OES D.L.'s are given as radial/axial view				

Thulium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 168.9342

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Tm}(\text{OH})_x(\text{H}_2\text{O})_y^{+3-x}$

69	168.934
1947	1.1
1545	
Tm	
[Xe]4f ¹³ 6s ²	2,3
9.32	

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, H₂SO₄, and HNO₃. Avoid HF, H₃PO₄, and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO₃/LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO₃/LDPE container.

Tm CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in H₂O/HNO₃); Ores (carbonate fusion in Pt⁰ followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 H₂O/HCl or HNO₃).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 313.126 nm	0.005/0.003 µg/mL	1	ion	U, Th, Be
ICP-OES 346.220 nm	0.008/0.006 µg/mL	1	ion	Rh, U
ICP-MS 169 amu	1 ppt	n/a	M+	¹⁵³ Eu ¹⁶ O
*ICP-OES D.L.'s are given as radial/axial view				

Ytterbium

LOCATION: Period 6 (lanthanoid)

ATOMIC WEIGHT: 173.054

COORDINATION NUMBER: 6 to 9, 10 for some compounds

CHEMICAL FORM IN SOLUTION: $\text{Yb}(\text{OH})_x(\text{H}_2\text{O})_y^{+3-x}$

70	173.054
1194	1.1
824	
Yb	
[Xe]4f ¹⁴ 6s ²	
6.97	2,3

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl , H_2SO_4 , and HNO_3 . Avoid HF , H_3PO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Yb CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in $\text{H}_2\text{O}/\text{HNO}_3$); Ores (carbonate fusion in Pt^0 followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 $\text{H}_2\text{O}/\text{HCl}$ or HNO_3).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 328.937 nm	0.002/0.0003 $\mu\text{g}/\text{mL}$	1	ion	U, Ce, V
ICP-OES 369.419 nm	0.003/0.0006 $\mu\text{g}/\text{mL}$	1	ion	Fe
ICP-MS 174 amu	2 ppt	n/a	M+	$^{158}\text{Gd}^{16}\text{O}$, $^{158}\text{Dy}^{16}\text{O}$, ^{174}Hf

*ICP-OES D.L.'s are given as radial/axial view

Lutetium

LOCATION: Group 13, Period 5

ATOMIC WEIGHT: 174.9668

COORDINATION NUMBER: 6

CHEMICAL FORM IN SOLUTION: $\text{Lu}(\text{OH})_x(\text{H}_2\text{O})_y^{+3-x}$

71	174.9668
3395	
1663	
Lu	
[Xe]4f ¹⁴ 5d ⁶ s ²	
9.84	1

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl, H₂SO₄, and HNO₃. Avoid HF, H₃PO₄, and neutral to basic media. Stable with most metals and inorganic anions forming an insoluble carbonate, oxide, oxalate, and fluoride. Avoid mixing with elements/solutions containing moderate amounts of fluoride.

STABILITY: 2–100 ppb levels stable for months in 1% HNO₃/LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO₃/LDPE container.

Lu CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in acids); Oxide (dissolved by heating in H₂O/HNO₃); Ores (carbonate fusion in Pt⁰ followed by HCl dissolution); Organic Matrices (dry ash and dissolve in 1:1 H₂O/HCl or HNO₃).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 261.542 nm	0.001/0.0003 µg/mL	1	ion	Th, Mo, V, W
ICP-OES 291.139 nm	0.006/0.0006 µg/mL	1	ion	Cr, U
ICP-MS 175 amu	1 ppt	n/a	M+	¹⁵⁹ Tb ¹⁶ O
*ICP-OES D.L.'s are given as radial/axial view				

Thorium

LOCATION: Period 7 (actinoid)

ATOMIC WEIGHT: 232.0381

COORDINATION NUMBER: 8

CHEMICAL FORM IN SOLUTION: $\text{Th}(\text{OH})^{3+}$ and $\text{Th}(\text{OH})_2^{2+}$

90
232.038

4788
1.1

1755

Th

[Rn]6d²7s²
4

11.7

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid H_3PO_4 , H_2SO_4 , and HF, although solubilities may not be a problem depending upon pH and matrix (i.e. ThF_4 is soluble in acids). Avoid neutral to basic media. Th^{4+} is stable with most metals and inorganic anions forming an insoluble carbonate, oxide, fluoride, oxalate, sulfate, and phosphate in neutral to slightly acidic media.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

Th CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (soluble in aqua regia); Oxide (the heated oxide is not soluble in acids except hot conc. H_2SO_4); Ores (Na_2O_2 fusion at $480 \pm 20^\circ\text{C}$ for 7 minutes, cool, and treat sintered mass with 50 mL cold water and let stand until disintegrated. The mass is transferred to a beaker and acidified with HCl, with 25 mL excess HCl added. Any residue is collected on a Whatman No. 42 filter, dried and ignited to 1000°C in Pt^0 crucible and ash treated with H_2SO_4 /HF and fumed. If residue remains, then treat it by peroxide fusion as described above).

ATOMIC SPECTROSCOPIC INFORMATION: (*italic text indicates severe at ~ concs.*)

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 283.730 nm	0.07/0.007 $\mu\text{g/mL}$	1	ion	U, Zr
ICP-OES 283.231 nm	0.07/0.007 $\mu\text{g/mL}$	1	ion	U, Mo, Ti, Fe, Cr
ICP-OES 274.716 nm	0.08/0.008 $\mu\text{g/mL}$	1	ion	Ti, Ta, <i>Fe</i> , V
ICP-MS 232 amu	1 ppt	n/a	M+	

*ICP-OES D.L.'s are given as radial/axial view

Uranium

LOCATION: Period 7 (actinoid)

ATOMIC WEIGHT: 238.0289

COORDINATION NUMBER: 8

CHEMICAL FORM IN SOLUTION: UO_2^{2+} (uranyl)

92
4134
1132
238.029
1.2

U

[Rn]5f³6d7s²
19.0
3,4,5,6

STORAGE & HANDLING: Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

CHEMICAL COMPATIBILITY: Soluble in HCl and HNO_3 . Avoid H_3PO_4 . H_2SO_4 and HF matrices should not be a problem depending upon [U]. Although the UO_2^{2+} ion is distinctly basic, any U^{+4} will ppt. in basic media. UO_2^{2+} salts are generally soluble in water and UO_2^{2+} is stable with most metals and inorganic anions. The uranyl phosphate is insoluble in water. UF_4 and UF_6 are water soluble.

STABILITY: 2–100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1–10,000 ppm solutions chemically stable for years in 2–5% HNO_3 /LDPE container.

U CONTAINING SAMPLES (PREPARATION & SOLUTION): Metal (dissolves rapidly in HCl and HNO_3); Oxide (soluble in HNO_3); Ores (digest for 1–2 hours with 1 gram of ore to 30 mL 1:1 HNO_3 . Silica insolubles are removed by filtration after bringing the sample to fumes with conc. H_2SO_4).

ATOMIC SPECTROSCOPIC INFORMATION:

TECHNIQUE / LINE	ESTIMATED D.L.*	ORDER	TYPE	INTERFERENCES
ICP-OES 385.958 nm	0.3/0.01 $\mu\text{g/mL}$	1	ion	Th, Fe
ICP-OES 367.007 nm	0.3/0.02 $\mu\text{g/mL}$	1	ion	Th, Ce
ICP-OES 263.553 nm	0.3/0.01 $\mu\text{g/mL}$	1	ion	Ce, Ir, Th, Rh, W, Zr, Ta, Ti, V, Hf, Fe, Re, Ru
ICP-MS 238 amu	2 ppt	n/a	M+	$^{206}\text{Pb}^{16}\text{O}_2$

*ICP-OES D.L.'s are given as radial/axial view



Dr. Paul R. Gaines and Madeline Gozzi discussing their research for the Osmium Primary Certified Reference Material (PCRM™).

Paul R. Gaines, Ph.D.

Dr. Paul R. Gaines has over four decades of spectroscopic experience. After earning his Ph.D. in chemistry at Iowa State University, Dr. Gaines worked in the laboratories of Exxon Research and Engineering and Union Carbide. Today, Dr. Gaines is the Senior Technical Advisor and Chairman of the Board of Directors at Inorganic Ventures, as well as an accomplished web author of many popular guides and papers for fellow spectroscopists.



Refine your results. Redefine your industry.

300 Technology Drive
Christiansburg, VA 24073 | USA

Phone: 1.540.585.3030

Fax: 1.540.585.3012

Toll-free: 1.800.669.6799

inorganicventures.com

✓ REGISTERED
ISO 9001
QSR-1034 BY QSR

✓ ACCREDITED
ISO 17025
0883.01 BY A2LA

✓ ACCREDITED
ISO 17034
0883.02 BY A2LA